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
## Agilent OpenLAB CDS 2.x od HPST, s.r.o.

Autor: Michal Novotný

Verze: 1.00 Leden 2018

---

# Obsah

---

<b>Obsah</b>	<b>2</b>
Úvodem . . . . .	3
OpenLAB Control Panel . . . . .	4
Instruments . . . . .	4
Projects . . . . .	4
Administration . . . . .	5
OpenLAB Acquisition . . . . .	6
Status . . . . .	6
Method . . . . .	6
Sequence . . . . .	7
Single Sample . . . . .	8
Práce s běžící analýzou . . . . .	8
OpenLAB Data Analysis . . . . .	8
Data Selection . . . . .	9
Data Processing . . . . .	9
Master metody vs. Sequence metody . . . . .	11
Práce s okny . . . . .	11
Integrace . . . . .	11
Identifikace . . . . .	13
Kalibrace . . . . .	14
System Suitability . . . . .	17
Specifická nastavení . . . . .	17
Single Injection report . . . . .	20
Sequence Summary Report . . . . .	21

---

# Úvodem

Milí zákazníci,

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
Roles	
Name	Description
Everything	All privileges
System Administrator	Manage users and security settings
Instrument Administrator	Manage instruments and locations
<b>Project Administrator</b>	<b>Manage projects and project groups</b>
Instrument User	View and run instruments
Technician	Laboratory technician
Chemist	Analytical chemist
Sample Scheduler for OpenLAB Administrator	Generic Role for Sample Scheduler Administrators
Sample Scheduler for OpenLAB Analyst	Generic Role for Sample Scheduler Analysts

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

Sequence run in progress – Prirucka123.rslt

Edit Exit View Review Completed Injections

General	State	Vial	Sample type	R
1	Acquiring	1	Blank	
2	Pending	2	Cal. Std.	Cle

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
Online Signals

Overlaid

- Grad. Pump: Pressure (bar)
- Column Oven: Temperature (°C)
- WWD: Absorbance (mAU) at 250 nm
- WWD: Signal without Reference (mAU) at 250 nm

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Sequence run in progress – Prirucka123.rslt

Edit Exit View Review Completed Injections

General	State	Vial	Sample type	Ru
1	Completed	1	Blank	
2	Completed	2	Cal. Std.	Clea
3	Acquiring	3	Cal. Std.	
4	Pending	4	Sample	

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



Processing Method window showing integration parameters for the DAD1A signal. The window is titled "Processing Method" and has a sidebar with categories: Properties, Signals, Extraction, Integration Events ChemStation, Compounds, and System Suitability. The "Integration Events ChemStation" section is expanded to show "Standard" parameters.

Global	Use	Time (min)	Event	Value
DAD1A	<input checked="" type="checkbox"/>	0.000	Slope sensitivity	1.000
	<input checked="" type="checkbox"/>	0.000	Peak width	0.020
	<input checked="" type="checkbox"/>	0.000	Area reject	1.000
	<input checked="" type="checkbox"/>	0.000	Height reject	1.700
	<input checked="" type="checkbox"/>	0.000	Shoulders mode	Off
	<input checked="" type="checkbox"/>	0.000	Area% reject	0.000

Buttons: Add integration event, Delete integration event

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
Chromatograms window with Chromatogram Properties dialog box open. The dialog box is titled "Chromatogram Properties" and has a sidebar with categories: Chromatograms, Integrator events, Spectrum extraction, Manual integration, Expected compounds, Compound windows, and Annotations. The "Integrator events" section is expanded to show "Change integrator timed events options".

**Change integrator timed events options**

- Display options
  - Show timed events
  - Restrict to focused injection
  - Show initial events
  - Show signal specific information
  - Show method information
  - Show annotation box
  - Annotation: Event name
- Tools
  - Edit timed events
  - Show command hints

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
$$Concentration = \frac{Amount * Multipliers * Dil. factors}{SampleAmount}$$

**Processing Method**

Properties  
Signals  
Extraction  
Chromatogram  
Spectrum  
Integration Events ChemStation  
Standard  
Advanced  
Manual Integration  
Compounds  
Identification  
**Calibration**  
Spectra

Compound Table **General**

External standard  Internal standard

Number of levels: 5

Curve calculation: From average per level

RF definition: Response per amount

**!** If you change the RF definition, you need to clear your calibration curve otherwise your results will be wrong. Use the "Clear all calibration" RunType for the 1st standard in the injection list to remove the old calibration curve.

Normalize to: 100.00 %  Include ISTD amount

Concentration calculation: Amount \* Multipliers \* Dil. factor  Calculate mass %

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**Processing Method**

Advanced  
Manual Integration  
Compounds  
Identification  
Calibration  
Spectra  
**System Suitability**  
Properties  
Column  
Tools  
Custom Calculation  
Reports

Column performance: All peaks

Pharmacopeias:  EP  USP  JP

Signal to noise: All peaks

Noise calculation: P2P

Noise range:  Automatic (On blank recommended for EP)  Fixed from: 2.00 to 3.00 min on current  Relative from: 0.00 to 0.40 min before peak start on blank

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### Processing Method

Advanced  
Manual Integration

- Compounds
- Identification
- Calibration
- Spectra
- System Suitability
- Properties
- Column
- Tools
- Custom Calculation
- Reports

Compound Table   **UV Impurity Check**   UV Confirmation

Calculate UV Purity    None  
 Identified peaks  
 **All peaks**

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### Processing Method

Advanced  
Manual Integration

- Compounds
- Identification
- Calibration
- Spectra
- System Suitability
- Properties
- Column
- Tools
- Custom Calculation
- Reports

Compound Table   UV Impurity Check   UV Confirmation

#	Type	Name	Impurity sensitivity [%]	UV ref
1	☺	Latka A	50	
2	☺	peak@4.225min	50	
3	☺	peak@5.056min	50	
4	☺	peak@5.893min	50	
5	☺	peak@6.687min	50	

Spectrum at RT: 4.071

94			
----	--	--	--


- Copy
- Paste
- Fill down
- Numeric precision...
- Add compound
- Add named group
- Add timed group
- Add/Update MS compound
- Delete selected compound
- Auto extract UV spectra
- Delete selected reference :

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
The screenshot shows the Agilent software interface with two main windows: 'Report Editor - Short\_Estd.rdl' and 'Report Preview'. The 'Report Editor' window displays a 'Single Injection Report' for 'Sample 1'. The report includes a header with the Agilent Technologies logo and a body with a table of metadata and a chromatogram. The metadata table is as follows:

Sample name:	Sample 1	Operator:	SYSTEM
Data file:	007_008-0301.D	Injection date:	2011-03-30 12:44:28+01:00
Instrument:	wat11166	Location:	Vial 8
Ini volume:	5.0	Tune:	Sample
Acq. method:	PERF_N.M	Calib Level:	
Processing method:	*Prirucka.nm	Sample amount:	30.00
Manually modified:	None		

The chromatogram shows a plot of mAU versus time with several peaks labeled: 'A1: 1st peak @ 2.53min', 'A2: 2nd peak @ 3.25min', 'A3: 3rd peak @ 3.83min', 'A4: 4th peak @ 4.38min', and 'A5: 5th peak @ 4.93min'. The 'Report Preview' window on the right shows a summary of the report data, including 'Sample name: Sample 1', 'Data file: 007\_008-0301.D', 'Instrument: wat11166', 'Inj. volume: 5.0', 'Acq. method: PERF\_N.M', and 'Processing method: \*Prirucka.nm'. It also shows a small chromatogram plot.

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
Header

Body

Sample name: Sample 1  
 Data file: 007\_008-0301.D  
 Instrument: wadt1166  
 Inj. volume: 5.0  
 Acq. method: PERF\_N.M  
 Processing method: \*Prirucka.pmx  
 Manually modified: None

Operator: SYSTEM  
 Injection date: 2011-03-30 12:44:28+01:00  
 Location: Vial 8  
 Type: Sample  
 Sample Multipliers: 1;1;1;1

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
**Chromatogram Report** 

Sample name: Sample 2  
 Data file: \_008\_008-0401.D  
 Instrument: wadt1166  
 Inj. volume: 6.0  
 Acq. method: PERF\_N.M  
 Processing method: \*Prirucka.pmx  
 Manually modified: None

Operator: SYSTEM  
 Injection date: 2011-03-30 12:58:08+01:00  
 Location: Vial 8  
 Type: Sample  
 Sample Multipliers: 1;1;1;1

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
The screenshot shows the 'Report Editor' window. The left sidebar contains 'Reporting' and 'Report Items' sections. The main area displays a 'System Suitability Report' for 'Krásná Fir'. The report includes a header with the company logo, a body with a table of instrument and sample information, a chromatogram plot, and a table of peak results.

**Instrument and Sample Information:**

Instrument:	wadt1166	Sample type:	Calibration
Acc. method:	PERF N.M	Acc. operator:	SYSTEM
Analysis method:	*Prirucka.omx		
Sequence Name	Example Data for ESTD		
Last changed:	2017-12-29 15:05:22+01:00		

**Chromatogram:** A plot showing detector response (mAU) versus time (min). The x-axis ranges from 0 to 12 minutes, and the y-axis ranges from 0 to 140 mAU. Several peaks are visible, with the most prominent ones between 4 and 7 minutes. The plot is titled 'DAD1 A, Sig=254,4 Ref=360,100'.

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Style

Advanced

### Customize the peaks and repeating properties

Repeat Table On

=Compound\_Name Repeat On Field Properties

Sort On:

=Compound\_Name Ascending

Set Filter On Peak Types

Identified Peaks

Compounds used for quantitation (Main Peak)

Missing Peaks

Unknown Peaks

Remove Duplicate Peak(s)

Area Rejection

Do not report peaks with area less than 0


Do not report peaks with area% less than 0

OK Cancel Apply

Compound Name	Latka A				
Sample Name	Name	RT	Area	Amount Unit	Concentration [%]
Calibration 3.2	Latka A	4.071	187.6883	34.457	34.4574
Calibration 3.1	Latka A	4.075	150.5901	30.196	30.1965
Compound Name	peak@4.225min				
Sample Name	Name	RT	Area	Amount Unit	Concentration [%]
Calibration 3.2	peak@4.2 25min	4.227	302.5599	34.491	34.4907
Calibration 3.1	peak@4.2 25min	4.231	242.2945	30.196	30.1956

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
Instrument: wadt1166      Sample type: Calibration  
 Acq. method: PERF\_N.M  
 Analysis method: \*Prirucka.pmx      Acq. operator: SYSTEM  
 Sequence Name: Example\_Data\_for\_ESTD  
 Last changed: 2017-12-29 15:05:22+01:00

Instrument: wadt1166      Sample type: Calibration  
 Acq. method: PERF\_N.M  
 Analysis method: \*Prirucka.pmx      Acq. operator: SYSTEM  
 Sequence Name: Example\_Data\_for\_ESTD  
 Last changed: 2017-12-29 15:05:22+01:00

Compound Name	Latka A				
Sample Name	Name	RT	Area	Amount Unit	Concentration [%]
Calibration 3.2	Latka A	4.071	187.6883	34.457	34.4574
Calibration 3.1	Latka A	4.075	150.5901	30.196	30.1965
		Average	169.1392		
		RSD	15.5094		

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## Calibration Table per Compounds

**Instrument:** wadt166 **Sample type:** Calibration  
**Acq. method:** PERF\_N.M **Acq. operator:** SYSTEM  
**Analysis method:** \*Prirucka.pmx  
**Sequence Name:** Example\_Data\_for\_ESTD  
**Last changed:** 2017-12-29 15:05:22+01:00

Compound Name	Latka A			
Sample Name	Level	Calib Amount	RT	Area
Calibration 1.1	1	10.000	4.082	37.3579
Calibration 1.2	1	10.000	4.068	49.8530
Calibration 2.1	2	20.000	4.069	86.8032
Calibration 2.2	2	20.000	4.068	112.2919
Calibration 3.1	3	30.000	4.075	150.5901
Calibration 3.2	3	30.000	4.071	187.6883


Compound Name	peak@4.225min			
Sample Name	Level	Calib Amount	RT	Area
Calibration 1.1	1	10.000	4.235	59.2781
Calibration 1.2	1	10.000	4.221	80.0910
Calibration 2.1	2	20.000	4.224	140.1595
Calibration 2.2	2	20.000	4.224	180.4457

## Calibration Table per Compounds

**Instrument:** wadt166 **Sample type:** Calibration  
**Acq. method:** PERF\_N.M **Acq. operator:** SYSTEM  
**Analysis method:** \*Prirucka.pmx  
**Sequence Name:** Example\_Data\_for\_ESTD  
**Last changed:** 2017-12-29 15:05:22+01:00

**Compound:** Latka A  
**Signal:** DAD1A  
**Exp. RT:** 4.069  
**Corr. Coeff.:** 0.994720  
**Residual:** 9.18985  
**RF RSD%:**  
**R<sup>2</sup>:** 0.98947  
**Formula:**  $y = ax + b$   
**a:** 5.63360  
**b:** -6.43090  
**c:** 0.00000  
**d:**

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