



# Agilent 1260 Infinity II Binary Pump

## Data Sheet



### Product Description

The Agilent 1260 Infinity II Binary Pump is ideally suited for a broad range of LC applications using conventional or superficially porous particles columns, even STM column technology can be used. Gradient formation is based on a high-pressure mixing principle and delay volumes (standard or low) can be easily configured by the user. The Agilent 1260 Infinity II Binary Pump is the pump of choice for reproducible gradients and high-performance providing high-throughput and fast separations including the use of STM column technology. For routine applications, with UV or MS detection, where high speed and resolution with uncompromised data quality are required, this is the pump of choice.

### Features

- Configurable delay volume - down to 120  $\mu\text{L}$  together with a flow range up to 5 mL/min provides universal applicability.
- Fast and easy - change from standard to low delay volume configuration is enabled.
- High gradient performance - even at low % B and narrow-bore flow rates.
- Integrated 2-channel-degasser.
- Fast and precise gradients - the pump is the perfect choice for fast and precise gradients using LC/MS, as well as UV-only systems.
- Fully exploits the speed and separation potential of Poroshell.



## Specifications

**Table 1** Physical Specifications

Type	Specification	Comments
Weight	17.6 kg (38.8 lbs)	
Dimensions (height × width × depth)	180 x 396 x 436 mm (7.1 x 15.6 x 17.2 inches)	
Line voltage	100 – 240 V~, ± 10 %	Wide-ranging capability
Line frequency	50 or 60 Hz, ± 5 %	
Power consumption	90 VA / 74 W	
Ambient operating temperature	4–55 °C (39–131 °F)	
Ambient non-operating temperature	-40 – 70 °C (-40 – 158 °F)	
Humidity	< 95 % r.h. at 40 °C (104 °F)	Non-condensing
Operating altitude	Up to 3000 m (9842 ft)	
Non-operating altitude	Up to 4600 m (15092 ft)	For storing the module
Safety standards: IEC, EN, CSA, UL	Installation category II, Pollution degree 2	For indoor use only.
ISM Classification	ISM Group 1 Class B	According to CISPR 11

**Table 2** Performance Specifications 1260 Infinity II Binary Pump (G7112B)

Type	Specification
Hydraulic system	Two dual piston in series pumps with servo-controlled variable stroke drive, power transmission by gears and ball screws, floating pistons
Flow range	settable: 0.001 – 5 mL/min recommended: 0.05 – 5.0 mL/min
Flow precision	≤0.07 % RSD or < 0.02 min SD, whichever is greater
Flow accuracy	± 1 % or 10 µL/min, whichever is greater
Pressure operating range	Up to 60 MPa (600 bar, 8702 psi) up to 5 mL/min
Pressure pulsation	< 2 % amplitude (typically < 1.3 %), or < 0.3 MPa (3 bar, 44 psi), whichever is greater <b>Low delay volume configuration:</b> < 5 % amplitude (typically < 2 %)
Compressibility compensation	Pre-defined, based on mobile phase compressibility
Recommended pH range	1.0 – 12.5
Gradient formation	High-pressure binary mixing
Delay volume	<b>Standard delay volume configuration:</b> 600 – 900 µL, (includes 400 µL mixer), dependent on back pressure <b>Low delay volume configuration:</b> 120 µL
Composition range	settable: 0 – 100 % recommended: 1 – 99 % or 5 µL/min per channel, whichever is greater
Composition precision	< 0.15 % RSD or < 0.04 min SD, whichever is greater
Composition accuracy	± 0.35 % absolute
Integrated degassing unit	Number of channels: 2 Internal volume per channel: 1.5 mL
Control	Agilent control software (e.g. ChemStation, EZChrom, OpenLAB CDS, MassHunter)
Instrument Control	Lab Advisor B.02.08 or above LC and CE Drivers A.02.14 or above
Local control	Agilent Instant Pilot (G4208A)
Communications	Controller-area network (CAN), Extended Remote Interface (ERI), Local Area Network (LAN)
Safety and maintenance	Extensive diagnostics, error detection and display through Agilent LabAdvisor, leak detection, safe leak handling, leak output signal for shutdown of the pumping system. Low voltage in major maintenance areas.

**Table 2** Performance Specifications 1260 Infinity II Binary Pump (G7112B)

Type	Specification
GLP features	Early maintenance feedback (EMF) for continuous tracking of instrument usage in terms of seal wear and volume of pumped mobile phase with pre-defined and user settable limits and feedback messages. Electronic records of maintenance and errors.
Housing	All materials are recyclable

## Ordering Details

**Table 3** Ordering Details 1260 Infinity II Binary Pump

Description	Product Number
<b>Agilent 1260 Infinity II Binary Pump</b> 600 bar pressure limit.	G7112B
HPLC System Tool Kit for 1260/1290 Infinity II LC	G7112B#001
Lab Advisor Advanced Software	G7112B#004
1260 Infinity II Max Uptime Kit	G7112B#007
1260 Infinity II Low Dispersion Kit MCT	G7112B#008
Active Seal Wash	G7112B#030
Solvent Selection Valve	G7112B#031
Delete Solvent Cabinet	G7112B#960

[www.agilent.com/chem/infinitylab-lc-series](http://www.agilent.com/chem/infinitylab-lc-series)

This information is subject to change without notice.

© Agilent Technologies, Inc., 2016  
Published November 1, 2016  
5991-7599EN