







# ClearSeq Comprehensive Cancer

## Confident Answers

The ClearSeq Comprehensive Cancer was designed for deep coverage of target bases and high on-target specificity, ensuring that variants are detected accurately with minimal sequencing (Figure 1). The research panel provides excellent performance even with challenging samples like formalin-fixed, paraffin-embedded (FFPE) DNA, which is important given that many solid tumor samples are available only as FFPE tissue (Figure 2). ClearSeq Comprehensive Cancer is compatible with the SureSelect Target Enrichment System and the combination of an optimized design and robust capture protocols allows the cancer researcher to confidently identify disease-associated somatic variants. The panel can also be easily customized by the addition of genes using SureDesign.




**Figure 1.** Excellent coverage of target bases when enriching HapMap cell line DNA with different SureSelect library preparation protocols (*XT* – post-capture library pooling; *QXT* – transposase-based library prep). Data is representative of eight replicates per protocol and each sample was sequenced to 300X depth (240 Mb).



**Figure 2.** Sequencing performance of SureSelect<sup>XT</sup> libraries from FFPE and fresh-frozen (FF) breast and lung tumor samples enriched using ClearSeq Comprehensive Cancer (240 Mb sequencing/sample). The level of DNA degradation is indicated by the DNA Integrity Number (DIN) provided by the 2200 TapeStation System, where a DIN of 10 and 1 indicate intact gDNA and completely degraded gDNA respectively.

## Comprehensive Solution

Sequence data generated from libraries enriched by ClearSeq Comprehensive Cancer are ideally suited for analysis using SureCall, an NGS data analysis software developed to address the critical bioinformatics needs of clinical researchers. Using guided cancer analysis workflows, identify known variants of interest within minutes, greatly reducing time-to-results. Together with robust sample QC and throughput scalability by automation, ClearSeq Comprehensive Cancer accelerates the profiling of many genes associated with cancer.



**Figure 3.** Accelerate time-to-results with a simplified 3-step analysis workflow using SureCall. Starting from FASTQ or BAM files, obtain a report of categorized variants within minutes.

**SureCall Mutation Report**

PREPARED BY: AGTDA02

SAMPLE ATTRIBUTES:  
Sample name: Sample\_42\_KnownVariant  
Date: 08 Jul 2014

Chr	Start	Stop	Allele	Type	HGVS(Genomic)
CHR1	43615072	43615073	A	SNP	NC_000001.10:g.43615072G>A
CHR1	115258747	115258748	T	SNP	NC_000001.10:g.115258747C>T
CHR9	5073770	5073771	T	SNP	NC_000009.11:g.5073770G>T
CHR17	74732785	74732786	T	SNP	17:74732785-T
CHR18	42531913	42531914	A	SNP	NC_000018.8:g.42531913G>A
CHR20	31023395	31023396	A	SNP	NC_000020.10:g.31023395G>A
CHR21	38164432	38164433	C	SNP	NC_000021.8:g.38164432T>C

DESCRIPTION/DETAILS: Found Variants

**NC\_000020.10:g.31023395G>A Category 1**

Pathology Lab Comment: Introduction of Stop Codon

Mutation Impact: This variant has effect on ASXL1, Transcript NM\_015338 of ASXL1 are affected. The Variant is HIGH impact and STOP\_GAINED for NM\_015338. This results in codon change of TgG/TgA and Amino acid change of W95D\*

## Ordering Information:

Part Number	Description
5190-8011	ClearSeq Comprehensive Cancer, 16, XT
5190-8012	ClearSeq Comprehensive Cancer, 96, XT
5190-8013	ClearSeq Comprehensive Cancer, 96 auto, XT
5190-8014	ClearSeq Comprehensive Cancer Plus, 16, XT
5190-8015	ClearSeq Comprehensive Cancer Plus, 96, XT
5190-8016	ClearSeq Comprehensive Cancer Plus, 96 auto, XT
5190-8017	ClearSeq Comprehensive Cancer, 16, XT2
5190-8018	ClearSeq Comprehensive Cancer, 96, XT2
5190-8019	ClearSeq Comprehensive Cancer, 96 auto, XT2
5190-8020	ClearSeq Comprehensive Cancer Plus, 16, XT2
5190-8021	ClearSeq Comprehensive Cancer Plus, 96, XT2
5190-8022	ClearSeq Comprehensive Cancer Plus, 96 auto, XT2



Find an Agilent customer center in your country:

[www.agilent.com/genomics/clearseqcancer](http://www.agilent.com/genomics/clearseqcancer)

U.S. and Canada: 1-800-227-9770 | [cag\\_sales-na@agilent.com](mailto:cag_sales-na@agilent.com)

For Research Use Only. Not for Use in Diagnostic Procedures.

© Agilent Technologies, Inc. 2015

Printed in USA, May 18, 2015

5991-5573EN



Agilent Technologies