

DNA Integrity Number (DIN) with the Agilent 2200 TapeStation System and the Agilent Genomic DNA ScreenTape Assay

Technical Overview

Authors

Marcus Gassmann
Agilent Technologies, Inc.
Waldbronn, Germany

Barry McHoull
Agilent Technologies, Inc.
Edinburgh, UK

Introduction

The results of numerous molecular screening and assay methods often rely on the overall quality of the genomic DNA (gDNA) input material. For example, array comparative genome hybridization (aCGH) and Next Generation Sequencing (NGS) can require intact, high quality gDNA to ensure high quality, unambiguous results. It is highly recommended to perform a quality control (QC) of the input material for these workflows, especially costly procedures such as NGS, as it saves both wasted time and effort on low quality samples.

The Agilent 2200 TapeStation system, in conjunction with the Agilent Genomic DNA ScreenTape assay, provides an excellent solution for assessing the quantity and integrity of the gDNA starting material. However, to provide a more direct measure of gDNA integrity, as well as to standardize integrity measurements, a software algorithm has been developed for the Genomic DNA ScreenTape assay. Based on the analysis of approximately 7,000 different gDNA samples, this algorithm provides a numerical assessment of gDNA integrity referred to as the DNA Integrity Number (DIN). This functionality, in conjunction with the TapeStation Analysis Software (revision A.01.05 or higher), automatically determines and displays DIN as a measure of gDNA integrity. Data files generated with a previous version of the software can be reanalyzed after a software upgrade to determine DIN.

This Technical Overview evaluates the robustness, in terms of reliability and reproducibility, of DIN generated with the 2200 TapeStation system and the Genomic DNA ScreenTape assay for gDNA samples, including Formalin-Fixed Paraffin-Embedded (FFPE) samples.



Agilent Technologies

Experimental

Materials

Commercially available intact human gDNA was purchased from Promega (Madison, WI, USA). A Digital Sonifier was used from Branson Ultrasonics (Danbury, CT, USA). Insulin syringes Omnican® 50 (30G × 8 mm) were obtained from B. Braun Melsungen AG (Melsungen, Germany). The 2200 TapeStation system (p/n G2965AA) with the Agilent TapeStation Analysis Software (revision A.01.05), Genomic DNA ScreenTape consumable (p/n 5067-5365) and Genomic DNA Reagents (p/n 5067-5366) were obtained from Agilent Technologies (Waldbronn, Germany).

Samples

Approximately 7,000 gDNA samples of eukaryotic origin including blood, dried blood spots, saliva, invertebrate, varying human tissues, and FFPE samples were kindly provided by customers. For robustness validation, commercially available gDNA was degraded using ultrasonication (10 % amplitudes, 2 second treatment, 5 second pause, up to 120 second treatment) or shearing with a fine gauge needle or a combination of both to generate a set of 15 gDNA samples with a wide range of gDNA degradation (gDNA Sample 1 highest integrity and gDNA Sample 15 highest degradation).

Genomic DNA analysis with the Agilent 2200 TapeStation system

DNA analysis was performed according to the Genomic DNA ScreenTape system Quick Guide¹. In brief, 1 µL gDNA was mixed with 10 µL Genomic DNA Sample buffer. Genomic DNA ScreenTape consumable, filtered loading tips, and the prepared samples were placed in the 2200 TapeStation instrument. The 2200 TapeStation system loaded, electrophoresed, imaged, and presented digitally analyzed results in less than 2 minutes per sample.

Results and Discussion

The gel image, as well as the electropherogram overlay in Figure 1, demonstrate the shift in DNA size and distribution with increasing degradation from right to left. Highly intact gDNA samples migrate as a single fragment in the gel image and a well-defined peak in the electropherogram above the largest ladder peak (48,500 bp). With increasing degradation, a shoulder of small gDNA fragments was formed and the main peak was shifted towards smaller sizes. Highly degraded gDNA appears as a smear in the gel image, and migrates as broad peak in the electropherogram with

sizes below 2,000 bp. Both gel images and electropherogram traces allow visual determination of a gDNA sample for integrity, however this approach is highly subjective. To provide an objective and standardized tool for reliable integrity assessment, the DNA Integrity Number (DIN) algorithm was developed. DIN determines the fragmentation of a genomic DNA sample by assessing the distribution of signal across the size range, and applies an automatically calculated number. DIN was generated using approximately 7,000 samples from a wide variety of genomic DNA sources comparing the signal against each and every sample.




Figure 1. A degradation series of 15 gDNA samples (at 60 ng/µL) was analyzed using the Agilent 2200 TapeStation system and the Agilent Genomic DNA ScreenTape assay to determine DIN. A) The gel image with DIN displayed below each lane for gDNA sample 15 to 1 from left to right. B) The electropherogram overlay of all 15 gDNA samples.

To provide a numerical assessment, the samples were assorted according to their signal distribution to a scale of DIN 1 to 10. A high DIN indicates highly intact gDNA, and a low DIN a strongly degraded gDNA sample.

Robustness and reproducibility

A set of 15 gDNA samples representative for the entire degradation range of the 7,000 samples was generated. The determined scores for the 15 gDNA samples ranged from DIN 9.8 to 1.3. Figure 1 provides a gel image and an overlay of the electropherogram for the 15 samples, and DIN is displayed directly under the individual lanes of the gel image. Using this standardized sample set the robustness and reproducibility of the algorithm in terms of independency of sample concentration, Genomic DNA Reagent batch, Genomic DNA ScreenTape device batch, TapeStation system type, and sample position on 96-well plates was analyzed.

Three genomic DNA samples representative for high, middle, and low integrity samples were analyzed at six different concentrations, covering the whole DIN specified functional range from 5 to 300 ng/μL, to demonstrate that DIN is independent of the sample concentration. Figure 2 clearly shows that the determined DIN for all three gDNA samples does not depend on the loaded DNA concentration within the DIN functional range.

A subsequent set of experiments was performed to prove that there is no influence of ScreenTape device batch (n = 4) or Genomic DNA Loading Buffer batch (n = 3) on the determined DIN for all 15 gDNA samples.

Two different instruments are available for the 2200 TapeStation system; the standard instrument is suitable for RNA, DNA, and protein analysis. The Nucleic Acid instrument is suitable for analysis of nucleic acids only. Both instruments types were used to analyze the set of 15 gDNA samples, to demonstrate that there is no effect of the instrument type on the determined DIN (data not shown).

The 2200 TapeStation system allows the analysis of samples from a 96-well plate. To test the stability and reproducibility of DIN, three gDNA samples selected across the integrity range were analyzed with n = 32 from a 96-well plate. Figure 3 clearly illustrates that the standard deviation for the average DIN for all three tested samples is small and, more importantly, stable across a 96-well plate measurement cycle.




Figure 2. Three gDNA samples were diluted to six different concentrations and analyzed using the Agilent 2200 TapeStation system and the Agilent Genomic DNA ScreenTape assay. Averaged DIN and standard deviation are plotted in this graph (n = 10 to 15).




Figure 3. Three gDNA samples (at 25 ng/μL) were analyzed from a 96-well plate using the Agilent 2200 TapeStation system and the Agilent Genomic DNA ScreenTape assay. The corresponding averaged DIN and the standard deviation are shown (n = 32).

Figure 4 summarizes the data from the above described experiments. The average DIN and the standard deviation were calculated and plotted.

DIN is highly reproducible, across different sample concentrations, instrument types, batches of ScreenTape device and Loading Buffer Reagent, and 96-well plate positions, respectively. The overall standard deviation was below 0.30 for all analyzed samples (Figure 4).

DIN availability for Genomic DNA TapeStation files generated with previous software versions

Taking advantage of the digital data obtained with the 2200 TapeStation system, the TapeStation Analysis Software automatically determines DIN (from revision A.01.05), and displays it underneath the gel image and in the sample data table. Data files generated with a previous software version A.01.05 can be reanalyzed after a software upgrade to determine DIN.

Conclusion

This Technical Overview demonstrates that DIN:

- Can measure a wide range of samples from highly intact gDNA (DIN 10) to very degraded gDNA (DIN ~1)
- Can be used to assess integrity of genomic DNA extracted from a variety of sample sources, including FFPE tissues
- Is robust throughout a wide concentration range
- Is highly reproducible across sample concentrations, instruments ScreenTape consumables and reagents, and 96-well plate positions




Figure 4. 15 gDNA samples were analyzed using the standard Agilent 2200 TapeStation system and the Agilent 2200 TapeStation Nucleic Acid system, using different ScreenTape and Agilent Genomic DNA Loading Buffer batches, as well as different sample concentrations. The graph and the table summarize the averaged DIN and standard deviation obtained combining the data of all performed measurements.

Reference

1. Agilent Genomic DNA ScreenTape System Quick Guide, *Agilent Technologies*, publication number G2964-90040 rev.C, 2014.

www.agilent.com/genomics/tapestation

This information is subject to change without notice.

© Agilent Technologies, Inc., 2014
Published in the USA, December 1, 2014
5991-5258EN