

Introduction

Cellular metabolism is a primary regulator of immune cell fate and function.

T cell activation is correlated with metabolic shifts, especially in glycolysis, providing increased metabolic flux to support high rates of cellular proliferation^{1,2}.

XF Hu T Cell Activation Assay Kit³ provides a rapid method to measure early activation-associated metabolic responses in human T cells, offering important insights into both the early dynamics of activation and its broader metabolic underpinning.

The assay employs a soluble CD3/CD28 activator, facilitating investigation of immune or metabolic modulators on T cell activation by providing a reliable and convenient workflow.

To exemplify the utility of the XF Hu T Cell Activation Assay in assessing the pharmacological modulation of activation, the effects of different tyrosine kinases on primary T cell activation were measured.

Experimental

XF Hu T Cell Activation Assay Workflow




Figure 1. Workflow summary XF Hu T Cell Activation Assay.

XF Hu T Cell Activation Assay Designs




Figure 2. Agilent Seahorse XF Hu T Cell Activation Assays measuring the proton efflux rate (PER) upon the injection of a soluble CD3/CD28 activator in the **standard assay** design (A) and the **modulation assay** design (B). The glycolytic rate can be analyzed by the absolute PER (pmol/min) or by the relative increase to the baseline, PER (%).

Data analysis strategy of XF Hu T Cell Activation Assay




Figure 3. Data analysis of Agilent Seahorse XF Hu T Cell Activation Assay using XF Seahorse Analytics. (A) and (B) are schematic illustrations for the methods calculating the **maximum rate** and the **AUC**, respectively.

Results and Discussion

XF Hu T Cell Activation Assay can rapidly measure T cell activation potential




Figure 4. Activation of naive CD4+ T cells by different amounts of the CD3/CD28 activator. The total injection volume of activator was fixed as 20 μL/well and the dose of activator was adjusted by diluting the indicated amount of activator in XF assay medium prior to loading in the injection port. (A) Assay kinetics exhibiting real-time changes in PER upon injection of the activator at different doses. (B) and (C) show maximum PER and AUC values at increasing activator doses. (D) IL-2 production at different activator doses as measured by ELISA with cells collected 72 hours post-activation.

Effective immune modulator dosing can be measured by using metrics calculated from T cell activation kinetic data




Figure 5. Modulator efficacy test using XF Hu T Cell Activation Assay Kit. Pre-activated CD4+ T cells were pretreated with different inhibitors at various concentrations for 30 min prior to the T Cell Activation Assay. (A) The kinetic graphs detected dose-dependent suppression of PER increase induced by CD3/CD28 Activator showing kinetic response variation depending on the inhibitor type. (B) The % maximum PER and % AUC values were calculated in Seahorse Analytics and plotted against the drug concentrations using GraphPad Prism to assess IC₅₀ values.

Acute effect of modulator can be evaluated by *in-situ* modulator injection




Figure 6. Evaluation of the acute effect of dasatinib and PD116326 on naive CD4+ T cell activation by using modulation assay protocol. Dasatinib and PD116326 were injected at the concentration indicated prior to the activator injection. The concentrations of dasatinib and PD116326 were determined based on the IC₅₀ value obtained from the data in Figure 5. (A) Kinetic graph of real-time changes in % PER; (B) % Maximum rates calculated; (C) % AUC calculated; (D) IL-2 production measured by ELISA using culture media collected at 72 hours post-cell activation.

Modulator effect on energy metabolism can be assessed by XF Real-time ATP Rate Assay⁴




Figure 7. Naive CD4+ T cells were pretreated with inhibitors and activated by beads conjugated with CD3/CD28 antibodies in Seahorse XF RPMI Medium, pH 7.4 containing 10 mM glucose, 1 mM pyruvate and 2 mM glutamine (200K cells/well). After 2 h of beads injection, ATP production rate was measured using the Seahorse XF Real-Time ATP Rate Assay (A). Accumulation of INF-γ (B) and IL-2 (C) were measured in the extracellular medium after 24 to 72 hours of activation.

Conclusions

By monitoring increases in PER, the XF platform enables a rapid method to visualize T cells activation in real-time.

The XF Hu T Cell Activation Assay Kit facilitates such measurements, allowing **early detection of T cell activation**, a real-time interrogation of activation kinetics, and an assessment of pharmacological or genetic activation modulation.

The assay employs a **soluble CD3/CD28 activator** and ready-to-use **PDL-coated XF Cell Culture Microplates** to maximize assay performance and provide a simple and convenient workflow.

Maximum PER provides an intuitive assessment of T cell activation potential, while AUC assessment enables interrogation of more complex kinetics.

Using the Standard assay design, these parameters can be used to generate quantitative **dose response metrics** useful for *in vitro* **assessments of immune modulators** as shown using Lck1 and Bcr-Abl kinase inhibitors.

Data is also presented illustrating how *in-situ* injections are utilized within the Modulation assay to conveniently reveal the immediate impact dasatinib in real-time, and these observations correlate with orthogonal ELISA measurements. This is particularly interesting as recent reports suggest dasatinib as a **potential CAR-T cell activation modulator**⁵.

Additional metabolic interrogations including **XF Real-time ATP Rate Assay** can also be incorporated to **determine the broader metabolic consequences of perturbed activation**.

These data therefore illustrate how the XF Hu T Cell Activation Assay can provide additional **insights into both chronic and acute immune modulators**, demonstrating its utility in cellular and molecular immunotherapy discovery for various diseases including cancer, immune dysfunction, and metabolic diseases.

Materials

- **T cells:** Naive CD4+ or CD8+ T cells (StemCell Technologies and Hemacare Corporation) and pre-activated cells cultured in ImmunoCult-XF T Cell Expansion Medium supplemented with IL-2 (StemCell Technologies) for 12 to 16 days after the activation were used.
- **Kits:** XF Hu T Cell Activation Kit (Agilent Technologies), XF Real-Time ATP Rate Assay Kit (Agilent Technologies).
- **Plate:** XF96 PDL Cell Culture Microplates (Agilent Technologies)
- **Modulators:**

Test Compound	Target	Test Compound	Target
Dasatinib	Bcr-Abl kinase	Imatinib	Bcr-Abl kinase
PD116326	Bcr-Abl kinase	TWS119	GSK-3β
Lck Inhibitor (CAS No. 213743-31-8)	Lck1 kinase	DCA	Pyruvate dehydrogenase

References

1. Gubser, P. M. *et al.* Rapid Effector Function of Memory CD8+ T Cells Requires an Immediate-Early Glycolytic Switch. *Nat. Immunol.* 2013, 14, 1064–1072.
2. Menk, A. V. *et al.* Early TCR Signaling Induces Rapid Aerobic Glycolysis Enabling Distinct Acute T Cell Effector Functions. *Cell Rep.* 2018, 22, 1509–1521.
3. Agilent Seahorse XF Hu T Cell Activation Assay Kit User Guide. Agilent Technologies technical overview, publication number: 5994-1811EN
4. Romero, N. *et al.* Quantifying Cellular ATP Production Rate Using Agilent Seahorse XF Technology. Agilent Technologies white paper, publication number: 5991-9303 EN.
5. Mestermann, K. *et al.* The Tyrosine Kinase Inhibitor Dasatinib Acts as a Pharmacological on/off Switch for CAR T Cells. *Sci. Transl. Med.* 2019, 11, eaau5907

For Research Use Only. Not for use in diagnostic procedures.