



QuikChange II-E Site-Directed Mutagenesis Kit

Instruction Manual

Catalog # 200555

Revision E.0

For Research Use Only. Not for use in diagnostic procedures.

200555-12



Agilent Technologies

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

Email

techservices@agilent.com

World Wide Web

www.genomics.agilent.com

Telephone

Location	Telephone
United States and Canada	800 227 9770
Austria	01 25125 6800
Benelux	02 404 92 22
Denmark	45 70 13 00 30
Finland	010 802 220
France	0810 446 446
Germany	0800 603 1000
Italy	800 012575
Netherlands	020 547 2600
Spain	901 11 68 90
Sweden	08 506 4 8960
Switzerland	0848 8035 60
UK/Ireland	0845 712 5292
All Other Countries	Please visit www.agilent.com/genomics/contactus

QuikChange II-E Site-Directed Mutagenesis Kit

CONTENTS

Materials Provided	1
Storage Conditions	1
Additional Materials Required	1
Notices To Purchaser	2
Introduction	3
QuikChange II-E Mutagenesis Control	5
Mutagenic Primer Design	6
Primer Design Guidelines.....	6
Additional Primer Considerations.....	7
Protocol	8
Mutant Strand Synthesis Reaction (Thermal Cycling).....	8
<i>Dpn</i> I Digestion of the Amplification Products.....	10
Transformation of XL1-Blue Electroporation-Competent Cells.....	10
Transformation Guidelines	13
Storage Conditions.....	13
Aliquoting Cells.....	13
Cuvette Gap Width.....	13
Preparing the Agar Plates for Color Screening.....	13
Troubleshooting	14
Preparation of Media and Reagents	15
References	15
MSDS Information	15
Quick-Reference Protocol	16

QuikChange II-E Site-Directed Mutagenesis Kit

MATERIALS PROVIDED

Materials provided	Quantity ^a
<i>PfuUltra</i> High-Fidelity DNA polymerase (2.5 U/ μ l)	25 U
10 \times reaction buffer	500 μ l
<i>Dpn</i> I restriction enzyme (10 U/ μ l)	100 U
Oligonucleotide control primer #1 [34-mer (100 ng/ μ l)] 5' CCA TGA TTA CGC CAA GCG CGC AAT TAA CCC TCA C 3'	750 ng
Oligonucleotide control primer #2 [34-mer (100 ng/ μ l)] 5' GTG AGG GTT AAT TGC GCG CTT GGC GTA ATC ATG G 3'	750 ng
pWhitescript 4.5-kb control plasmid (5 ng/ μ l)	50 ng
dNTP mix ^{b,c}	10 μ l
XL1-Blue electroporation-competent cells ^d (blue tubes)	5 \times 100 μ l
pUC18 control plasmid (0.1 ng/ μ l in TE buffer ^e)	10 μ l
StrataClean resin	100 μ l

^a Reagents are provided for 10 total reactions, which includes 5 control reactions.

^b Thaw the dNTP mix once, prepare single-use aliquots, and store the aliquots at -20°C . **Do not subject the dNTP mix to multiple freeze-thaw cycles.**

^c **The composition of the dNTP mix is proprietary.** This reagent has been optimized for the QuikChange II-E site-directed mutagenesis protocols and has been qualified for use in conjunction with the other kit components. Do not substitute with dNTP mixes provided with other Agilent kits.

^d Genotype: *recA1 endA1 gyrA96 thi-1 hsdR17 supE44 relA1 lac* [F' *proAB lac^aZ Δ M15 Tn10* (Tet^r)]

^e See *Preparation of Media and Reagents*.

STORAGE CONDITIONS

XL1-Blue Electroporation-Competent Cells and pUC18 Control Plasmid: -80°C

All Other Components: -20°C

ADDITIONAL MATERIALS REQUIRED

5-Bromo-4-chloro-3-indolyl- β -D-galactopyranoside (X-gal)

Isopropyl-1-thio- β -D-galactopyranoside (IPTG)

14-ml BD Falcon polypropylene round-bottom tubes (BD Biosciences catalog #352059)

NOTICES TO PURCHASER

Use of this product is licensed under one or more of the following U.S. Patent Nos. 5,789,166, 5,932,419, 6,391,548, 6,713,285, 7,132,265, and 7,176,004.

Notice to Purchaser: Limited License

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,258,569, 6,171,785, 6,127,155, 6,030,787, 5,994,056, 5,876,930, 5,804,375, 5,789,224, 5,773,258 (claims 1 and 6 only), 5,723,591, 5,677,152 (claims 1 to 23 only), 5,618,711, 5,538,848, and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA

INTRODUCTION

In vitro site-directed mutagenesis is an invaluable technique for characterizing the dynamic, complex relationships between protein structure and function, for studying gene expression elements, and for carrying out vector modification. Several approaches to this technique have been published, but these methods generally require single-stranded DNA (ssDNA) as the template¹⁻⁴ and are labor intensive or technically difficult. Our QuikChange II-E Site-Directed Mutagenesis Kit* allows site-specific mutation in virtually any double-stranded plasmid, thus eliminating the need for subcloning and for ssDNA rescue.⁵ In addition, the QuikChange II-E site-directed mutagenesis kit does not require specialized vectors, unique restriction sites, multiple transformations or *in vitro* methylation treatment steps. The rapid three-step procedure generates mutants with greater than 80% efficiency in a single reaction (see Figure 1). The protocol is simple and uses either miniprep plasmid DNA or cesium-chloride-purified DNA. For long (~8 kb) or difficult targets, we offer the QuikChange II XL Site Directed Mutagenesis Kits (Catalog #200521 and #200522).




FIGURE 1 Overview of the QuikChange II-E site-directed mutagenesis method.

* U.S. Patent Nos. 6,391,548, 5,932,419, 5,789,166, 7,132,265, and 7,176,004.

The QuikChange II-E site-directed mutagenesis kit is used to make point mutations, replace amino acids, and delete or insert single or multiple adjacent amino acids. The QuikChange II-E site-directed mutagenesis method is performed using *PfuUltra* high-fidelity (HF) DNA polymerase** for mutagenic primer-directed replication of both plasmid strands with the highest fidelity.^{||} The basic procedure utilizes a supercoiled double-stranded DNA (dsDNA) vector with an insert of interest and two synthetic oligonucleotide primers, both containing the desired mutation (see Figure 1). The oligonucleotide primers, each complementary to opposite strands of the vector, are extended during temperature cycling by *PfuUltra* HF DNA polymerase, without primer displacement. Extension of the oligonucleotide primers generates a mutated plasmid containing staggered nicks. Following temperature cycling, the product is treated with *Dpn* I. The *Dpn* I endonuclease (target sequence: 5'-Gm⁶ATC-3') is specific for methylated and hemimethylated DNA and is used to digest the parental DNA template and to select for mutation-containing synthesized DNA.⁶ (DNA isolated from almost all *E. coli* strains is *dam* methylated and therefore susceptible to *Dpn* I digestion.) The nicked vector DNA containing the desired mutations is purified with StrataClean resin and then transformed into XL1-Blue electroporation-competent cells.

Note *While plasmid DNA isolated from almost all of the commonly used E. coli strains (dam⁺) is methylated and is a suitable template for mutagenesis, plasmid DNA isolated from the exceptional dam⁻ E. coli strains, including JM110 and SCS110, is not suitable.*

Unwanted second-site errors are virtually eliminated and high mutation efficiencies are obtained using this method due to the high fidelity of the *PfuUltra* HF DNA polymerase, the use of a small amount of starting DNA template and the use of a low number of thermal cycles.

** U.S. Patent Nos. 6,734,293, 6,444,428, 6,183,997, and 5,948,663.

^{||} *PfuUltra* HF DNA polymerase has 18-fold higher fidelity in DNA synthesis than *Taq* DNA polymerase.

QUIKCHANGE II-E MUTAGENESIS CONTROL

The pWhitescript 4.5-kb control plasmid is used to test the efficiency of mutant plasmid generation using the QuikChange II-E site-directed mutagenesis kit. The pWhitescript 4.5-kb control plasmid contains a stop codon (TAA) at the position where a glutamine codon (CAA) would normally appear in the β -galactosidase gene of the pBluescript II SK(-) phagemid (corresponding to amino acid 9 of the protein). XL1-Blue electroporation-competent cells transformed with this control plasmid appear white on LB-ampicillin agar plates (see *Preparation of Media and Reagents*), containing IPTG and X-gal, because β -galactosidase activity has been obliterated. The oligonucleotide control primers create a point mutation on the pWhitescript 4.5-kb control plasmid that reverts the T residue of the stop codon (TAA) at amino acid 9 of the β -galactosidase gene to a C residue, to produce the glutamine codon (CAA) found in the wild-type sequence. Following transformation, colonies can be screened for the β -galactosidase (β -gal⁺) phenotype of blue color on media containing IPTG and X-gal.

MUTAGENIC PRIMER DESIGN

Note *Mutagenic primers can be designed using our web-based QuikChange Primer Design Program available online at www.agilent.com/genomics/qcpd.*

Primer Design Guidelines

The mutagenic oligonucleotide primers for use in this protocol must be designed individually according to the desired mutation. The following considerations should be made when designing mutagenic primers:

- ◆ Both of the mutagenic primers must contain the desired mutation and anneal to the same sequence on opposite strands of the plasmid.
- ◆ Primers should be between 25 and 45 bases in length, with a melting temperature (T_m) of $\geq 78^\circ\text{C}$. Primers longer than 45 bases may be used, but using longer primers increases the likelihood of secondary structure formation, which may affect the efficiency of the mutagenesis reaction.
- ◆ The following formula is commonly used for estimating the T_m of primers:

$$T_m = 81.5 + 0.41(\% \text{GC}) - (675/N) - \% \text{ mismatch}$$

For calculating T_m :

- N is the primer length in bases
- values for **%GC** and **% mismatch** are whole numbers

For calculating T_m for primers intended to introduce insertions or deletions, use this modified version of the above formula:

$$T_m = 81.5 + 0.41(\% \text{GC}) - (675/N)$$

where N does not include the bases which are being inserted or deleted.

Note *When using primer design software for QuikChange site-directed mutagenesis applications, be aware that the T_m calculated by the primer design software may differ from the T_m value calculated using the formula presented above. We recommend verifying primer T_m 's using the formula above or by using the QuikChange T_m calculator, available online at www.genomics.agilent.com.*

- ◆ The desired mutation (deletion or insertion) should be in the middle of the primer with ~10–15 bases of correct sequence on both sides.
- ◆ The primers optimally should have a minimum GC content of 40% and should terminate in one or more C or G bases.

Additional Primer Considerations

- ◆ The mutagenesis protocol uses 125 ng of each oligonucleotide primer. To convert nanograms to picomoles of oligo, use the following equation:

$$X \text{ pmoles of oligo} = \frac{\text{ng of oligo}}{330 \times \# \text{ of bases in oligo}} \times 1000$$

For example, for 125 ng of a 25-mer:

$$\frac{125 \text{ ng of oligo}}{330 \times 25 \text{ bases}} \times 1000 = 15 \text{ pmole}$$

- ◆ Primers need not be 5' phosphorylated but **must** be purified either by fast polynucleotide liquid chromatography (FPLC) or by polyacrylamide gel electrophoresis (PAGE). Failure to purify the primers results in a significant decrease in mutation efficiency.
- ◆ It is important to keep primer concentration in excess. We suggest varying the amount of template while keeping the concentration of the primer constantly in excess.

PROTOCOL

Mutant Strand Synthesis Reaction (Thermal Cycling)

Notes *Ensure that the plasmid DNA template is isolated from a dam⁺ E. coli strain. The majority of the commonly used E. coli strains are dam⁺. Plasmid DNA isolated from dam⁻ strains (e.g. JM110 and SCS110) is not suitable.*

To maximize temperature-cycling performance, we strongly recommend using thin-walled tubes, which ensure ideal contact with the temperature cycler's heat blocks. The following protocols were optimized using thin-walled tubes.

1. Synthesize two complimentary oligonucleotides containing the desired mutation, flanked by unmodified nucleotide sequence. Purify these oligonucleotide primers prior to use in the following steps (see *Mutagenic Primer Design*).

2. Prepare the control reaction as indicated below:

- 5 μ l of 10 \times reaction buffer
- 2 μ l (10 ng) of pWhitescript 4.5-kb control plasmid (5 ng/ μ l)
- 1.25 μ l (125 ng) of oligonucleotide control primer #1
[34-mer (100 ng/ μ l)]
- 1.25 μ l (125 ng) of oligonucleotide control primer #2
[34-mer (100 ng/ μ l)]
- 1 μ l of dNTP mix
- 38.5 μ l ddH₂O (to bring the final reaction volume to 50 μ l)

Then add

- 1 μ l of *PfuUltra* HF DNA polymerase (2.5 U/ μ l)

3. Prepare the sample reaction(s) as indicated below:

Note *Set up a series of sample reactions using various amounts of dsDNA template ranging from 5 to 50 ng (e.g., 5, 10, 20, and 50 ng of dsDNA template) while keeping the primer concentration constant.*

- 5 μ l of 10 \times reaction buffer
- X μ l (5–50 ng) of dsDNA template
- X μ l (125 ng) of oligonucleotide primer #1
- X μ l (125 ng) of oligonucleotide primer #2
- 1 μ l of dNTP mix
- ddH₂O to a final volume of 50 μ l

Then add

- 1 μ l of *PfuUltra* HF DNA polymerase (2.5 U/ μ l)

- Cycle each reaction using the cycling parameters outlined in Table I. (For the control reaction, use a 5-minute extension time and run the reaction for 12 cycles.)

TABLE I

Cycling Parameters for the QuikChange II-E Site-Directed Mutagenesis Method

Segment	Cycles	Temperature	Time
1	1	95°C	30 seconds
2	12–18	95°C	30 seconds
		55°C	1 minute
		68°C	1 minute/kb of plasmid length*

* For example, a 5-kb plasmid requires 5 minutes at 68°C per cycle.

- Adjust segment 2 of the cycling parameters according to the type of mutation desired (see the following table):

Type of mutation desired	Number of cycles
Point mutations	12
Single amino acid changes	16
Multiple amino acid deletions or insertions	18

- Following temperature cycling, place the reaction on ice for 2 minutes to cool the reaction to $\leq 37^{\circ}\text{C}$.

Note *If desired, amplification may be checked by electrophoresis of 10 μl of the product on a 1% agarose gel. A band may or may not be visualized at this stage. In either case proceed with Dpn I digestion and transformation.*

Dpn I Digestion of the Amplification Products

1. Add 1 μl of the Dpn I restriction enzyme (10 U/ μl) directly to each amplification reaction.
2. Gently and thoroughly mix each reaction mixture by pipetting the solution up and down several times. Spin down the reaction mixtures in a microcentrifuge for 1 minute and immediately incubate each reaction at 37°C for 1 hour to digest the parental (i.e., the nonmutated) supercoiled dsDNA.

Transformation of XL1-Blue Electroporation-Competent Cells

Notes *Please read the Transformation Guidelines before proceeding with the transformation protocol.*

XL1-Blue cells are resistant to tetracycline. If the mutagenized plasmid contains only the tet^R resistance marker, an alternative tetracycline-sensitive strain of competent cells must be used.

Pre-Electroporation Clean-Up using StrataClean Resin

1. Vortex the StrataClean resin to completely resuspend it.
2. Using a sterile pipet tip, pipet 7 μl of the slurry into the tube containing the mutagenesis reaction.
3. Vortex this mixture for 15 seconds, then incubate the mixture at room temperature for 1 minute.
4. Pellet the resin by centrifugation at $>2000 \times g$ for 1 minute.
5. Using a sterile pipet tip, carefully remove the supernatant containing the DNA to a fresh microcentrifuge tube, taking care to avoid the resin pellet.

Note *It is sufficient to collect most (~80%) of the supernatant volume, since only a small volume is used in the transformation reaction. DO NOT attempt to remove all of the supernatant, as this may result in carry-over of resin that could interfere with subsequent procedures.*

Electroporation of XL1-Blue Competent Cells

1. Before beginning, ensure that the electroporation cuvettes (0.1-cm gap), and 1.5-ml microcentrifuge tubes have been thoroughly chilled on ice. Preheat sterile SOC medium (see *Preparation of Media and Reagents*) to 37°C.

2. Set the electroporator to the following settings:

Applied volts (0.1-cm gap)	Field strength (0.1-cm gap)	Resistance	Capacitance
1700 V	17 kV/cm	200 Ω	25 μ F

3. Thaw the electroporation-competent cells on ice.
4. Aliquot 40 μ l of cells into each of three chilled 1.5-ml microcentrifuge tubes (one for the experimental mutagenesis, one for the control mutagenesis, and one for the pUC18 transformation control).

Note *Electroporation-competent cells are very fragile. Pipet gently when aliquoting the cells.*

5. For the experimental and control mutagenesis reactions, add 2 μ l of resin-purified mutagenesis reaction to 40 μ l of cells. Gently mix the cells and DNA and keep the mixture on ice.

Dilute the pUC18 control plasmid 1:10 with sterile, distilled water (dH₂O). Add 1 μ l of the diluted pUC18 control plasmid to 40 μ l of cells. Gently mix the cells and DNA and keep the mixture on ice.

6. Transfer each DNA-cell mixture to a **chilled** electroporation cuvette (0.1-cm gap), tapping the cuvette until the mixture settles evenly to the bottom.
7. Slide the cuvette into the chilled electroporation chamber until the cuvette connects with the electrical contacts.
8. Pulse the sample once, then quickly remove the cuvette. **Immediately** add 960 μ l of sterile SOC medium (held at 37°C) to resuspend the cells.
9. Transfer each cell suspension to a sterile 14-ml BD Falcon polypropylene round-bottom tube. Incubate the tubes at 37°C for 1 hour while shaking at 225–250 rpm.

10. Plate the appropriate volume of each transformation reaction, as indicated in the table below.

Plate the experimental mutagenesis reaction on agar plates containing the appropriate antibiotic for the plasmid vector.

For the mutagenesis and transformation controls, spread cells on LB–ampicillin agar plates containing 80 µg/ml X-gal and 20 mM IPTG (see *Preparing the Agar Plates for Color Screening in Transformation Guidelines*).

Transformation Reaction Plating Volumes

Reaction Type	Volume to Plate
pWhitescript mutagenesis control	100 µl
pUC18 transformation control	2.5 µl (in 100 µl of SOC)*
Sample mutagenesis	100 µl on each of two plates

* Place a 100-µl pool of SOC broth on the agar plate, pipet the 2.5 µl of the transformation reaction into the pool, then spread the mixture.

11. Incubate the transformation plates at 37°C for >16 hours.

Expected Results for the Control Transformations

The expected colony number from the transformation of the pWhitescript control mutagenesis reaction is ≥50 colonies. Greater than 80% of the colonies should contain the mutation and appear as blue colonies on agar plates containing IPTG and X-gal.

Note The mutagenesis efficiency (ME) for the pWhitescript 4.5-kb control plasmid is calculated by the following formula:

$$ME = \frac{\text{Number of blue colony forming units (cfu)}}{\text{Total number of colony forming units (cfu)}} \times 100\%$$

If transformation of the pUC18 control plasmid was performed, >250 colonies should be observed (transformation efficiency >10¹⁰ cfu/µg pUC18) with >98% of the colonies having the blue phenotype.

Expected Results for Sample Transformations

The expected colony number is between 10 and 500 colonies, depending upon the base composition and length of the DNA template employed. For suggestions on increasing colony number, see *Troubleshooting*. The insert of interest should be sequenced to verify that selected clones contain the desired mutation(s).

TRANSFORMATION GUIDELINES

Storage Conditions

The XL1-Blue electroporation-competent cells are very sensitive to even small variations in temperature and must be stored at the bottom of a -80°C freezer. Transferring tubes from one freezer to another may result in a loss of efficiency. The cells should be placed at -80°C directly from the dry ice shipping container.

Aliquoting Cells

When aliquoting, keep the XL1-Blue electroporation-competent cells on ice at all times. It is essential that the microcentrifuge tubes and electroporation cuvettes are placed on ice before the cells are thawed and that the cells are aliquoted directly into the prechilled tubes.

Cuvette Gap Width

Use a cuvette with a 0.1-cm gap to maximize the transformation efficiency and to minimize the possibility of arcing. A cuvette with a 0.2-cm gap is not recommended because the transformation efficiency is lower and the possibility of arcing is higher.

Preparing the Agar Plates for Color Screening

To prepare the LB agar plates for blue–white color screening, add 80 $\mu\text{g}/\text{ml}$ of 5-bromo-4-chloro-3-indolyl- β -D-galactopyranoside (X-gal), 20 mM isopropyl-1-thio- β -D-galactopyranoside (IPTG), and the appropriate antibiotic to the LB agar. Alternatively, 100 μl of 10 mM IPTG and 100 μl of 2% X-gal can be spread on the LB agar plates 30 minutes prior to plating the transformations. Prepare the IPTG in sterile dH_2O ; prepare the X-gal in dimethylformamide (DMF). Do not mix the IPTG and the X-gal before pipetting them onto the plates because these chemicals may precipitate.

TROUBLESHOOTING

When used according to the guidelines outlined in this instruction manual, this kit provides a reliable means to conduct site-directed mutagenesis using dsDNA templates. Variations in the base composition and length of the DNA template and in thermal cycler performance may contribute to differences in mutagenesis efficiency. We provide the following guidelines for troubleshooting these variations

Observation	Suggestion(s)
Low transformation efficiency or low colony number	Ensure that sufficient mutant DNA is synthesized in the reaction. Increase the amount of the <i>Dpn</i> I-treated, resin-purified DNA used in the transformation reaction to 4 μ l.
	Visualize the DNA template on a gel to verify the quantity and quality. Nicked or linearized plasmid DNA will not generate complete circular product. Verify that the template DNA is at least 80% supercoiled.
	It is not uncommon to observe low numbers of colonies, especially when generating large mutations. Most of the colonies that do appear, however, will contain mutagenized plasmid.
	Ethanol precipitate the <i>Dpn</i> I digested PCR product, and resuspend in a decreased volume of water before transformation.
Low mutagenesis efficiency or low colony number with the control reaction	Different thermal cyclers may contribute to variations in ramping efficiencies. Adjust the cycling parameters for the control reaction and repeat the protocol for the sample reactions.
	Ensure that the electroporation-competent cells are stored at the bottom of a -80°C freezer immediately upon arrival. (see also <i>Transformation Guidelines</i>).
	Verify that the agar plates were prepared correctly. See <i>Preparing the Agar Plates for Color Screening</i> , and follow the recommendations for IPTG and X-Gal concentrations carefully.
	For best visualization of the blue (β -gal ⁺) phenotype, the control plates must be incubated for at least 16 hours at 37°C .
	Avoid multiple freeze-thaw cycles for the dNTP mix. Thaw the dNTP mix once, prepare single-use aliquots, and store the aliquots at -20°C . Do not subject the dNTP mix to multiple freeze-thaw cycles.
Low mutagenesis efficiency with the sample reaction(s)	Allow sufficient time for the <i>Dpn</i> I to completely digest the parental template; repeat the digestion if too much DNA template was present.
	Avoid multiple freeze-thaw cycles for the dNTP mix. Thaw the dNTP mix once, prepare single-use aliquots, and store the aliquots at -20°C . Do not subject the dNTP mix to multiple freeze-thaw cycles.
	The formation of secondary structures may be inhibiting the mutagenesis reaction. Increasing the annealing temperature up to 68°C may help to alleviate secondary structure formation and improve mutagenesis efficiency.
False positives	Poor quality primers can lead to false positives. Radiolabel the primers and check for degradation on an acrylamide gel or resynthesize the primers.
	False priming can lead to false positives. Increase the stringency of the reaction by increasing the annealing temperature up to 68°C .
Unwanted deletion or recombination of plasmid DNA following mutagenesis and transformation	Transform the mutagenesis reaction into competent cells that are designed to prevent recombination events, such as Agilent's SURE 2 Supercompetent Cells (Catalog #200152). Note that SURE 2 competent cells are not recommended for use with mutagenized plasmids greater than 10 kb in size; note also that SURE 2 cells are Kan ^r , Tet ^r , and Chl ^r , and are not compatible with plasmid selection using kanamycin, tetracycline, or chloramphenicol resistance markers.

PREPARATION OF MEDIA AND REAGENTS

<p>LB Agar (per Liter) 10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Add deionized H₂O to a final volume of 1 liter Adjust pH to 7.0 with 5 N NaOH Autoclave Pour into petri dishes (~25 ml/100-mm plate)</p>	<p>LB–Ampicillin Agar (per Liter) (Use for reduced satellite colony formation) 1 liter of LB agar Autoclave Cool to 55°C Add 100 mg of filter-sterilized ampicillin Pour into petri dishes (~25 ml/100-mm plate)</p>
<p>SOB Medium (per Liter) 20.0 g of tryptone 5.0 g of yeast extract 0.5 g of NaCl Add deionized H₂O to a final volume of 1 liter Autoclave Add 10 ml of filter-sterilized 1 M MgCl₂ and 10 ml of filter-sterilized 1 M MgSO₄ prior to use</p>	<p>SOC Medium (per 100 ml)</p> <p>Note <i>This medium should be prepared immediately before use.</i></p> <p>2 ml of filter-sterilized 20% (w/v) glucose or 1 ml of filter-sterilized 2 M glucose SOB medium (autoclaved) to a final volume of 100 ml</p>
	<p>TE Buffer 10 mM Tris-HCl (pH 7.5) 1 mM EDTA</p>

REFERENCES

1. Kunkel, T. A. (1985) *Proc Natl Acad Sci U S A* 82(2):488-92.
2. Sugimoto, M., Esaki, N., Tanaka, H. and Soda, K. (1989) *Anal Biochem* 179(2):309-11.
3. Taylor, J. W., Ott, J. and Eckstein, F. (1985) *Nucleic Acids Res* 13(24):8765-85.
4. Vandeyar, M. A., Weiner, M. P., Hutton, C. J. and Batt, C. A. (1988) *Gene* 65(1):129-33.
5. Papworth, C., Bauer, J. C., Braman, J. and Wright, D. A. (1996) *Strategies* 9(3):3-4.
6. Nelson, M. and McClelland, M. (1992) *Methods Enzymol* 216:279-303.

MSDS INFORMATION

Material Safety Data Sheets (MSDSs) are provided online at <http://www.genomics.agilent.com>. MSDS documents are not included with product shipments.

QuikChange II-E Site-Directed Mutagenesis Kit

Catalog #200555

QUICK-REFERENCE PROTOCOL

- ◆ Prepare the control and sample reaction(s) as indicated below:

Note Set up a series of sample reactions using various amounts of dsDNA template (e.g., 5, 10, 20, and 50 ng of dsDNA template).

Control Reaction

5 μ l of 10 \times reaction buffer
2 μ l (10 ng) of pWhitescript 4.5-kb control template (5 ng/ μ l)
1.25 μ l (125 ng) of oligonucleotide control primer #1 [34-mer (100 ng/ μ l)]
1.25 μ l (125 ng) of oligonucleotide control primer #2 [34-mer (100 ng/ μ l)]
1 μ l of dNTP mix
38.5 μ l ddH₂O (for a final volume of 50 μ l)

Sample Reaction

5 μ l of 10 \times reaction buffer
X μ l (5–50 ng) of dsDNA template
X μ l (125 ng) of oligonucleotide primer #1
X μ l (125 ng) of oligonucleotide primer #2
1 μ l of dNTP mix
ddH₂O to a final volume of 50 μ l

- ◆ Add 1 μ l of *PfuUltra* HF DNA polymerase (2.5 U/ μ l) to each control and sample reaction.
- ◆ Cycle the reactions using the cycling parameters outlined in the following table:

Segment	Cycles	Temperature	Time
1	1	95°C	30 seconds
2	12–18	95°C	30 seconds
		55°C	1 minute
		68°C	1 minute/kb of plasmid length

- ◆ Adjust segment 2 of the cycling parameters in accordance with the type of mutation desired (see the table in step 6 of *Mutant Strand Synthesis Reaction (Thermal Cycling)* in the instruction manual).
- ◆ Add 1 μ l of the *Dpn* I restriction enzyme (10 U/ μ l).
- ◆ Gently and thoroughly mix each reaction, spin down in a microcentrifuge for 1 minute, and immediately incubate at 37°C for 1 hour to digest the parental supercoiled dsDNA.

- ◆ Resuspend the StrataClean resin by vortexing, then add 7 μl of the suspension to each 50- μl reaction. Vortex the reactions for 15 seconds, incubate at room temperature for 1 minute, and then spin the reactions for 1 minute at $>2000 \times g$ to collect the resin.
- ◆ Transfer most of the supernatant to a fresh tube, taking care to avoid the resin.
- ◆ Transform 2 μl of the resin-purified DNA from each control and sample reaction into separate 40- μl aliquots of XL1-Blue electroporation-competent cells (see *Electroporation of XL1-Blue Competent Cells* in the instruction manual for a more detailed electroporation protocol).
- ◆ After outgrowth for 1 hour at 37°C, spread 100- μl samples of the transformed cell suspension onto LB agar plates containing the appropriate antibiotic for selection of the mutagenized plasmid. For the pWhitescript control, spread 100 μl of cells on LB-ampicillin agar plates containing 80 $\mu\text{g}/\text{ml}$ X-gal and 20 mM IPTG.
- ◆ Incubate the transformation plates at 37°C for >16 hours.