



# ***TaqPlus Precision PCR System***

## **Instruction Manual**

**Catalog #600210 (100 U), #600211 (500 U), and #600212 (1000 U)**

Revision D.0

**For Research Use Only. Not for use in diagnostic procedures.**

600210-12



**Agilent Technologies**

## LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

## ORDERING INFORMATION AND TECHNICAL SERVICES

### Email

techservices@agilent.com

### World Wide Web

www.genomics.agilent.com

### Telephone

Location	Telephone
United States and Canada	800 227 9770
Austria	01 25125 6800
Benelux	02 404 92 22
Denmark	45 70 13 00 30
Finland	010 802 220
France	0810 446 446
Germany	0800 603 1000
Italy	800 012575
Netherlands	020 547 2600
Spain	901 11 68 90
Sweden	08 506 4 8960
Switzerland	0848 8035 60
UK/Ireland	0845 712 5292
All Other Countries	Please visit <a href="http://www.agilent.com/genomics/contactus">www.agilent.com/genomics/contactus</a>

# **TaqPlus Precision PCR System**

## **CONTENTS**

<b>Materials Provided</b> .....	<b>1</b>
<b>Storage Conditions</b> .....	<b>1</b>
<b>Additional Materials Required</b> .....	<b>1</b>
<b>Notices to Purchaser</b> .....	<b>1</b>
<b>Introduction</b> .....	<b>2</b>
<b>Protocol Considerations</b> .....	<b>3</b>
<b>Polymerase Chain Reaction Protocol Using the <i>TaqPlus Precision PCR System</i></b> .....	<b>4</b>
<b>Troubleshooting</b> .....	<b>6</b>
<b>References</b> .....	<b>7</b>
<b>MSDS Information</b> .....	<b>7</b>

# TaqPlus Precision PCR System

## MATERIALS PROVIDED

Materials provided	Quantity		
	Catalog #600210	Catalog #600211	Catalog #600212
TaqPlus Precision polymerase mixture (5 U/ $\mu$ l)	100 U	500 U	1000 U
TaqPlus Precision 10 $\times$ buffer <sup>a</sup>	1 ml	2 $\times$ 1 ml	4 $\times$ 1 ml

<sup>a</sup> The TaqPlus Precision 10 $\times$  buffer provided with this system is optimized to produce PCR products with a minimal number of errors. Modifying the TaqPlus Precision 10 $\times$  buffer may adversely affect the fidelity and yield of the PCR product.

## STORAGE CONDITIONS

All Components:  $-20^{\circ}\text{C}$

## ADDITIONAL MATERIALS REQUIRED

Temperature cycler  
PCR tubes  
PCR primers  
Deoxynucleotides

## NOTICES TO PURCHASER

Limited Label License: The purchase of this product conveys to the purchaser only the limited, non-transferable right to use the purchased quantity of the product for the purchaser's own research by the purchaser only under patent claims outside the U.S. corresponding to the expired U.S. Patent No. 5,436,149 (claims 6-16). No rights are granted to the purchaser to sell, modify for resale or otherwise transfer this product, either alone or as a component of another product, to any third party. Takara Bio reserves all other rights, and this product may not be used in any manner other than as provided herein. For information on obtaining a license to use this product for purposes other than research, please contact Takara Bio Inc., Seta 3-4-1, Otsu, Shiga 520-2193, Japan (Fax +81-77-543-9254).

### NOTICE TO PURCHASER: LIMITED LICENSE

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,258,569, 6,171,785, 6,127,155, 6,030,787, 5,994,056, 5,876,930, 5,804,375, 5,789,224, 5,773,258 (claims 1 and 6 only), 5,723,591, 5,677,152 (claims 1 to 23 only), 5,618,711, 5,538,848, and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

## INTRODUCTION

---

The *TaqPlus Precision* PCR system<sup>1</sup> consists of the *TaqPlus Precision* polymerase mixture and the *TaqPlus Precision* 10× buffer, which are specifically formulated for high-fidelity polymerase chain reaction (PCR) applications. The *TaqPlus Precision* polymerase mixture is an optimized blend of Agilent's cloned *Pfu* DNA polymerase<sup>2</sup> and *Taq2000* DNA polymerase, a highly purified, recombinant version of *Taq* DNA polymerase. The *TaqPlus Precision* PCR system exhibits the highest replication accuracy achieved using DNA polymerase mixtures.<sup>3</sup>

Thermostable DNA polymerases offer distinct advantages for specific PCR applications. *Pfu* DNA polymerase exhibits the lowest error rate of any thermostable DNA polymerase analyzed.<sup>1, 4-6</sup> This proofreading enzyme is also characterized by a relatively slow polymerization rate requiring the use of long extension times (i.e.,  $\geq 2$  minutes/kb of DNA target) to generate high yields of PCR product.<sup>7</sup> *Taq* DNA polymerase exhibits a fivefold higher polymerization rate than *Pfu* DNA polymerase and is used routinely to generate high yields of PCR product  $\leq 5$  kb in length using short extension times (i.e.,  $\sim 1$  minute/kb of DNA target). *Taq* DNA polymerase possesses a significantly higher error rate than *Pfu* DNA polymerase. However, Agilent's *Taq2000* DNA polymerase is superior to other commercial *Taq* DNA polymerases in that this recombinant version of *Taq* DNA polymerase minimizes artifactual smearing in PCR amplification reactions requiring long extension times.<sup>8</sup>

DNA polymerase mixtures, consisting of a nonproofreading DNA polymerase and a proofreading DNA polymerase, are typically employed to amplify longer DNA targets and to generate higher yields of PCR product than can be obtained using either DNA polymerase alone.<sup>9-11</sup> Agilent offers two mixtures of *Taq* and *Pfu* DNA polymerases designed for PCR: the *TaqPlus Long* PCR system and the *TaqPlus Precision* PCR system. The *TaqPlus Long* PCR system typically generates the highest yields of PCR product and can be used to amplify unusually difficult or long targets (up to 35 kb). The *TaqPlus Precision* PCR system also generates relatively high yields of PCR product using short PCR extension times. The error rate, however, is significantly lower than the *TaqPlus Long* PCR system. In addition, the *TaqPlus Precision* PCR system successfully amplifies plasmid and lambda DNA templates up to 15 kb in length and single-copy genomic DNA templates up to 10 kb in length as well as templates that are difficult to amplify using single enzyme formulations.<sup>1</sup>

## PROTOCOL CONSIDERATIONS

---

- ◆ The *TaqPlus Precision* PCR system successfully amplifies single-copy genomic DNA templates up to 10 kb in length and plasmid and lambda DNA templates up to 15 kb in length.<sup>1</sup> Amplification of certain longer targets, especially single-copy chromosomal templates  $\geq 5$  kb, may require additional optimization of the reaction conditions (see *Troubleshooting*).
- ◆ As with all PCR amplification reactions, the following factors are critical for success: template quality, primer design, denaturation conditions, annealing temperature, and extension times.
- ◆ If using existing primer pairs, the annealing temperatures may require optimization in the *TaqPlus Precision* 10× buffer provided. If new primers are being designed, design the primers to have similar melting temperatures (typically between 60° and 80°C).
- ◆ The *TaqPlus Long* PCR system typically amplifies longer targets with higher yields than can be obtained with the *TaqPlus Precision* PCR system. If high fidelity is not a primary concern, we recommend the use of the *TaqPlus Long* PCR system to generate robust PCR products up to 35 kb in length.
- ◆ Cycling parameters require optimization specific for the primer–template sets and thermal cycler used. These optimized cycling parameters are not necessarily transferable between thermal cyclers designed by different manufacturers.

# POLYMERASE CHAIN REACTION PROTOCOL USING THE TAQPLUS PRECISION PCR SYSTEM

---

1. Prepare a reaction mixture on ice for the appropriate number of samples to be amplified. Add the components *in order* and mix gently. Table I provides an example of a reaction mixture for the amplification of a typical single-copy chromosomal target. The recipe listed in Table I is for one reaction and must be adjusted for multiple samples. The final volume of each sample reaction is 50  $\mu$ l.
2. Immediately before thermal cycling, aliquot 50  $\mu$ l of the reaction mixture into the appropriate number of sterile thin-wall PCR tubes or standard 0.5-ml microcentrifuge tubes and place the tubes on ice.

**TABLE I**

**Reaction Mixture for a Typical Single-Copy Chromosomal Locus PCR Amplification**

Component	Amount per reaction
Distilled water (dH <sub>2</sub> O)	41.1 $\mu$ l
<i>TaqPlus Precision</i> 10 $\times$ buffer <sup>a,b</sup>	5.0 $\mu$ l
100 mM dNTP mix (25 mM of each dNTP)	0.4 $\mu$ l
DNA template (250 ng/ $\mu$ l) <sup>c</sup>	1.0 $\mu$ l <sup>c</sup>
Primer #1 (100 ng/ $\mu$ l) <sup>d,e</sup>	1.0 $\mu$ l
Primer #2 (100 ng/ $\mu$ l) <sup>d,e</sup>	1.0 $\mu$ l
<i>TaqPlus Precision</i> polymerase mixture (5 U/ $\mu$ l) <sup>f</sup>	0.5 $\mu$ l (2.5 U)
Total reaction volume	50 $\mu$ l

- <sup>a</sup> The *TaqPlus Precision* 10 $\times$  buffer provided with this system is optimized to produce PCR products with a minimal number of errors. Modifying the reaction buffer may adversely affect the fidelity and yield of the PCR product.
- <sup>b</sup> If the *TaqPlus Precision* 10 $\times$  buffer appears cloudy, heat the buffer for ~5 minutes at 42°C until clear.
- <sup>c</sup> The amount of DNA template required will vary depending on the type of DNA being amplified. Generally 50–1000 ng of genomic DNA template is recommended; however, less DNA template can be used for amplification of lambda or plasmid PCR targets or for amplification of multicopy chromosomal genes. Use of higher amounts of DNA template reduces the number of PCR cycles required and therefore reduces the percentage of PCR products exhibiting mutations.
- <sup>d</sup> Primer concentrations between 0.3 and 0.5  $\mu$ M are recommended (generally 100–250 ng for typical 20- to 30-mer oligonucleotide primers in a 50- $\mu$ l reaction volume).
- <sup>e</sup> Design primers to have similar melting temperatures (typically between 60°C and 80°C).
- <sup>f</sup> The amount of *TaqPlus Precision* polymerase mixture will vary depending on the complexity and length of the template to be amplified. A typical starting point is 2.5 U per 50- $\mu$ l reaction. The standard range for most plasmid and lambda DNA templates is 0.5–1.0  $\mu$ l (2.5–5.0 U) and for most single-copy genomic DNA templates is 0.5–1.5  $\mu$ l (2.5–7.5 U) (see *Troubleshooting*).

**TABLE II**  
**PCR Cycling Parameters<sup>a,b</sup>**

Segment	Number of cycles	Temperature	Duration
1	1	95°C	2 minutes
2	30–35	95°C	40 seconds
		Primer $T_m - 5^\circ\text{C}$ <sup>c</sup>	30 seconds
		72°C	1 minute/kb PCR target
3	1	72°C	10 minutes

<sup>a</sup> Thin-wall PCR tubes are highly recommended. These PCR tubes are optimized to ensure more efficient heat transfer and to maximize thermal-cycling performance.

<sup>b</sup> Optimized cycling parameters are not necessarily transferable between thermal cyclers designed by different manufacturers; therefore, each manufacturer's recommendations for optimal cycling parameters should be consulted.

<sup>c</sup> The annealing temperature may require optimization. Typically annealing temperatures will range between 55° and 72°C.

3. If the temperature cycler is not equipped with a heated cover, overlay each reaction with ~50  $\mu\text{l}$  of DNase-, RNase-, and protease-free mineral oil (Sigma, St. Louis, Missouri).
4. Perform PCR using optimized cycling conditions. Suggested cycling parameters for the *TaqPlus Precision* PCR system are indicated in Table II.
5. Analyze the PCR amplification products by electrophoresis on a 1.0–4.0% (w/v) agarose gel.



## TROUBLESHOOTING

Observations	Suggestions
No product or low yield	Increase the amount of <i>TaqPlus Precision</i> polymerase mixture. This is especially important for plasmid and lambda targets >10 kb in length or for single-copy, genomic DNA targets >5 kb in length.
	Increase the amount of full-length intact DNA template and/or increase the number of cycles up to a maximum of 40 cycles.
	Use intact and highly purified DNA template.
	Use minimum primer concentrations of 0.3 $\mu\text{M}$ (i.e., at least 100 ng for typical 20-mer oligonucleotide primers in a 50- $\mu\text{l}$ reaction volume).
	Lower the annealing temperature in 5°C increments.
	Allow at least 60 seconds of extension time for each kilobase to be amplified (90 seconds of extension time per kilobase may also be helpful for difficult templates).
	Denaturation times of 30–60 seconds at 94–95°C are usually sufficient, while longer denaturation times may damage the DNA template; use the shortest denaturation time compatible with successful PCR on the thermal cycler.
	Primer pairs with matched primer melting temperatures ( $T_m$ ) and complete complementarity between the primer and template are recommended.
	Gel-purified primers $\geq 18$ nucleotides in length are desired for successful PCR.
	Adjust the ratio of primer versus template to optimize yield of the desired product.
Artifactual smears	Decrease the amount of <i>TaqPlus Precision</i> polymerase mixture.
	Reduce the extension time.
	Optimize the cycling parameters specifically for the primer–template set and the thermal cycler used.
Multiple bands	Increase the annealing temperature in 5°C increments and/or use a hot start. <sup>12, 13</sup>
	Use Perfect Match PCR enhancer to improve PCR product specificity or increase the annealing temperature as indicated previously.

## REFERENCES

---

1. Nielson, K. B., Cline, J., Bai, F., McMullan, D., McGowan, B. *et al.* (1997) *Strategies* 10(1):29–32.
2. Cline, J., Nielson, K. B., Scott, B. and Mathur, E. (1992) *Strategies* 5(2):50.
3. Stratagene (1996). A Guide to Pfu DNA Polymerase. Stratagene Cloning Systems, La Jolla, California.
4. Flaman, J. M., Frebourg, T., Moreau, V., Charbonnier, F., Martin, C. *et al.* (1994) *Nucleic Acids Res* 22(15):3259-60.
5. Cha, R. S. and Tilly, W. G. (1995). In *PCR Primer*, C. W. Dieffenbach and G. S. Dveksler (Eds.). Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York.
6. Cline, J., Braman, J. C. and Hogrefe, H. H. (1996) *Nucleic Acids Res* 24(18):3546-51.
7. Nielson, K. B., Costa, G. L. and Braman, J. (1996) *Strategies* 9(1):24–25.
8. Nielson, K. B., Scott, B., McMullan, D., Fowler, R., McGowan, B. *et al.* (1996) *Strategies* 9(1):7–8.
9. Barnes, W. (1994) *Proc Natl Acad Sci U S A* 91:2216–2220.
10. Cohen, J. (1994) *Science* 263:1564–1565.
11. Cheng, S., *et al.* (1994) *Proc Natl Acad Sci U S A* 91:5695–5699.
12. Chou, Q., Russell, M., Birch, D. E., Raymond, J. and Bloch, W. (1992) *Nucleic Acids Res* 20(7):1717-23.
13. Sharkey, D. J., Scalice, E. R., Christy, K. G., Jr., Atwood, S. M. and Daiss, J. L. (1994) *Biotechnology (N Y)* 12(5):506-9.

## MSDS INFORMATION

---

Material Safety Data Sheets (MSDSs) are provided online at <http://www.genomics.agilent.com>. MSDS documents are not included with product shipments.