


Seeding Adherent Cells in Agilent Seahorse XF96 Cell Culture Microplates

Basic Procedure

Agilent Seahorse XF Assays are performed in an Agilent Seahorse 96-well XF Cell Culture Microplate in conjunction with an XF°96/XF96 Sensor Cartridge. Each microplate is formatted in a typical 96-well design, as shown. The seeding surface of each well is 0.106 cm², smaller than in a



typical 96-well plate, but larger than in a typical 384 well plate. This procedure describes recommendations for seeding adherent cell types for use with the Agilent Seahorse XFe96/XF96 Analyzer.

- Harvest and re-suspend the cells to desired final concentration to seed in 80 μL of growth medium. Optimal cell seeding numbers vary widely, though are typically between 5,000 – 40,000 cells per well and must be determined empirically. For further information on optimal cell density, please see the Agilent Seahorse Cell Reference Database (http:// www.agilent.com/cell-reference-database/) and/or XF Assay Guides and Templates http://www.agilent.com/en-us/support/cell-analysis-(seahorse)/seahorse-assay-guides-templates.
- 2. Seed 80 µL of cell suspension per well (as shown in figure below); do not seed cells in background correction wells (A1, A12, H1, H12). Be sure to put medium only (no cells) in the background correction wells.
- 3. Optional: Allow plate to rest at room temperature in the tissue culture hood for one hour. This can promote even cell distribution and reduce edge effects for some cell types¹. Monitor adherence using a microscope.
- Allow the cells to grow overnight in a cell culture incubator. Monitor growth and health of cells using a microscope.

Hint: Hold the pipette tip at an angle about halfway down the side of the wells for best technique and most homogeneous cell layer.

¹. Lundholt BK, Scudder KM, Pagliaro L. A simple technique for reducing edge effect in cell-based assays. J Biomol Screen. 2003 Oct;8(5):566-7

Learn more

www.aqilent.com/en-us/promotions/seahorse-xf-technology

For Research Use Only. Not for use in diagnostic procedures.

Product specifications and descriptions in this document are subject to change without notice.

© Agilent Technologies, Inc., 2014, 2017 Published in the USA, June 30, 2017

