



MBS Mammalian Transfection Kit

Instruction Manual

Catalog #200388

Revision C.0

For Research Use Only. Not for use in diagnostic procedures.

200388-12



Agilent Technologies

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

Email

techservices@agilent.com

World Wide Web

www.genomics.agilent.com

Telephone

Location	Telephone
United States and Canada	800 227 9770
Austria	01 25125 6800
Benelux	02 404 92 22
Denmark	45 70 13 00 30
Finland	010 802 220
France	0810 446 446
Germany	0800 603 1000
Italy	800 012575
Netherlands	020 547 2600
Spain	901 11 68 90
Sweden	08 506 4 8960
Switzerland	0848 8035 60
UK/Ireland	0845 712 5292
All Other Countries	Please visit www.agilent.com/genomics/contactus

MBS Mammalian Transfection Kit

CONTENTS

Materials provided	1
Additional Materials Required	1
Storage Conditions	1
Introduction	2
96-Well Microtiter Plate Format	3
Transfection Protocol	3
Preparing the DNA Suspension for Transfection	4
Preparing the Cells for Transfection	4
Applying the DNA Suspension	4
Incubating the Transfection Reactions	5
Microtransfection Protocol Using 96-Well Microtiter Plates	5
Preparing the DNA Suspension for Microtransfection	5
Preparing the Cells for Microtransfection	6
Applying the DNA Suspension and Incubating the Transfection Reactions	6
β-Galactosidase Assay for 60-mm Dishes	7
Harvesting the Transfected Cells	7
Performing the β -Galactosidase Assay	7
β-Galactosidase Assay for 96-Well Microtiter Plates	8
β-Galactosidase Histochemical Staining Assay	9
General Tissue Culture Techniques: Transfection of Mammalian Cells	10
Culture Inoculation	10
Troubleshooting	10
Preparation of Media and Reagents	11
References	12
MSDS Information	12
Quick-Reference Protocols	13
Transfection Protocol	13
Microtransfection Protocol	13

MBS Transfection Kit

MATERIALS PROVIDED

Material provided	Composition	Quantity
Solution 1	2.5 M CaCl ₂	10 ml
Solution 2	2× BBS ^o (pH 6.95)	50 ml
Solution 3	Modified bovine serum (MBS)	80 ml
Control plasmid	pCMV β-gal	40 μg

^o BBS, *N,N*-bis(2-hydroxyethyl)-2-aminoethanesulfonic acid and buffered saline.

ADDITIONAL MATERIALS REQUIRED

Phosphate-buffered saline (PBS) (see *Preparation of Media and Reagents*)
Potassium ferricyanide
Potassium ferrocyanide
Formaldehyde and glutaraldehyde
5-Bromo-4-chloro-3-indolyl-β-D-galactopyranoside (X-gal)
5-ml BD Falcon polystyrene round bottom tubes (BD Biosciences catalog #352054)

STORAGE CONDITIONS

All Components: -20°C

NOTICE TO PURCHASER

The use of the CMV Promoter is covered under U.S. Patent Nos. 5,168,062 and 5,385,839 owned by the University of Iowa Research Foundation and licensed FOR RESEARCH USE ONLY. For further information, please contact UIRF at 319-335-4546.

Revision C.0

© Agilent Technologies, Inc. 2015.

INTRODUCTION

The CaPO₄ method of DNA transfection was initially developed as a technique for measuring the infectivity of isolated viral DNA.¹ The precipitate formed by CaPO₄ enables transfection by enhancing the adsorption of DNA to cell membranes, thus facilitating the DNA uptake by mammalian cells. The CaPO₄ precipitate also limits the digestion of DNA by DNase associated with mammalian cells.

The CaPO₄ technique has been adapted to introduce DNA from a variety of sources into mammalian cells. Initial attempts to produce stably transfected clones relied on the conversion of mutant cell lines, deficient in a selectable trait [e.g., thymidine kinase (TK)] to a wild-type phenotype.²⁻⁴ Later, DNA sequences not associated with vital functions were stably introduced into mammalian cells by either physically linking them to, or cotransfecting them with, the TK gene.^{5, 6} The utility of the CaPO₄ precipitation method was greatly enhanced by the development of the neomycin antibiotic selection technique, which allows for the selection of transfected cells by the acquisition of resistance to the neomycin analog, G418.

Transfection has been used for the isolation of functional genes and the analysis of transcriptional regulatory elements and of translational and RNA-processing signals. Most traditional techniques produce a low frequency of transfection, ranging from 1 in 10⁶ to 1 in 10⁴ cells.⁷⁻¹² The Agilent MBS Mammalian Transfection Kit*,¹³ is based on a technique that transiently transfects cells at a frequency of greater than 1 in 10 cells.¹⁴

All the necessary reagents and protocols for efficient, rapid and simplified DNA transfection of mammalian cells are included. The kit has been optimized for 100 transfections in 60-mm tissue culture dishes using a modified CaPO₄ protocol for the production of transient transfectants. Solutions 1 and 2 are used in the CaPO₄ precipitation protocol. Solution 3 is a specially modified bovine serum (MBS) solution that is added to the culture medium instead of fetal calf serum (FCS) during the incubation of the cells with the CaPO₄-DNA precipitate. This method is particularly suited for the generation of transient transfectants due to the short incubation time required; however, stable transfectants have also been generated with increased efficiency.

* U.S. Patent No. 5,589,368.

96-Well Microtiter Plate Format

The MBS Mammalian Transfection Kit can also be used for transfections performed using a 96-well microtiter plate format (see *Microtransfection Protocol Using 96-Well Microtiter Plates*). This protocol includes modifications to accommodate a large number of transfections carried out simultaneously in microtiter plates. Due to the high transfection efficiency, only small numbers of cells per well are required to obtain a significant signal from a transfected reporter gene such as β -galactosidase.

This protocol variation not only decreases the quantity of cells, DNA, transfection reagents and plasticware required, but also reduces the time required to carry out transient transfection assays. Cells can be lysed directly in the 96-well microtiter plate and the reporter gene product can be assayed in the plate, thus significantly reducing the manipulations required to carry out a large number of transfection assays. Applications of this method for performing multiple microtransfections include promoter analysis of numerous constructs or genetic screening protocols that aim at identifying agents which modulate gene expression.

TRANSFECTION PROTOCOL

Note *Please read the entire transfection or microtransfection protocol before beginning. See the final page of this instruction manual for condensed, quick-reference protocols that can be detached and used when detailed instructions are no longer required.*

This protocol outlines the procedure necessary for transient transfection using 60-mm tissue culture dishes. If using 100-mm tissue culture dishes, double the quantities and the volumes of the solutions included in the protocol.

1. Inoculate 60-mm tissue culture dishes with $1-5 \times 10^5$ exponentially growing cells per dish 24 hours before the transfection.
2. Grow the cells overnight in 5 ml of the appropriate culture medium.

Important *Roswell Park Memorial Institute (RPMI) media cannot be used for CaPO_4 transfection. The excess positive charge in this media will cause formation of a dense precipitate that is toxic to the cells.*

Preparing the DNA Suspension for Transfection

The optimal amount of DNA must be determined for each plasmid and cell line separately. The recommended amount is between 1 and 10 µg of DNA per dish. Start with 5 µg of control DNA per dish when preparing the pCMV β-gal control plasmid.

The amounts indicated in the following steps are for duplicate 60-mm tissue culture dishes (500 µl each).

1. For each duplicate transfection, dilute the optimal amount of DNA (2–20 µg) with sterile water to a final volume of 450 µl in a 5-ml BD Falcon polystyrene round bottom tube.

Note *Polystyrene tubes yield superior results compared to tubes composed of polyethylene or polypropylene.*

2. Add 50 µl of solution 1 and 500 µl of solution 2. Mix the contents of the tube by gently tapping the sides.
3. Incubate the DNA suspension at room temperature for 10–20 minutes. During this incubation period, prepare the cells as indicated in the following section (see *Preparing the Cells for Transfection*).

Preparing the Cells for Transfection

1. Remove the standard culture medium from the tissue culture dishes by aspiration.
2. Wash the cells twice with phosphate-buffered saline (PBS).
3. Add 5 ml of culture medium containing 4–7% solution 3 (i.e., MBS).

Note *The optimal concentration of MBS should be determined by titration. A good starting point for most cell lines is 6%.*

The cells are now ready to receive the DNA.

Applying the DNA Suspension

1. **Gently** resuspend any precipitate in the DNA suspension by pipetting the suspension up and down with a pipettor set at 500 µl. The DNA suspension should appear clear to opaque.
2. Slowly add 500 µl of the DNA suspension to each tissue culture dish **dropwise in a circular motion** to distribute the DNA suspension evenly onto each tissue culture dish.
3. Swirl each tissue culture dish **once**.

Incubating the Transfection Reactions

1. Incubate the dishes for 3 hours at 37°C under 5% CO₂.

Note *Incubation at 37°C under 5% CO₂ works well; however, for optimal transfection efficiency, the 3-hour incubation should be performed at 35°C under 3% CO₂.*

2. After 3 hours, check the precipitate, which will vary in consistency from slightly perceptible to noticeably granular.
3. Remove the culture medium by aspiration and wash the tissue culture dishes three times with PBS, which should be added gently to avoid dislodging the cells. Check for residual precipitate by microscopic inspection. Add additional wash steps as required to ensure complete removal of the precipitate.
4. Add 5 ml of standard culture medium and incubate the dishes overnight at 37°C under 5% CO₂.
5. Harvest the cells. For transient transfections, perform an assay to detect the gene product.

For the pCMV β-gal control plasmid, perform an assay to detect the β-galactosidase gene product in order to monitor the transfection efficiency. See either the *β-Galactosidase Assay for 60-mm Dishes* section or the *β-Galactosidase Histochemical Staining Assay* section.

MICROTRANSFECTION PROTOCOL USING 96-WELL MICROTITER PLATES

1. Inoculate a 96-well microtiter plate with 1–2 × 10⁴ cells per well in 0.2 ml of standard culture medium 24 hours before the transfection. The optimal cell density per well should be determined by seeding various amounts of cells.

Preparing the DNA Suspension for Microtransfection

1. Dilute 5–20 μg of DNA with sterile water to a final volume of 450 μl in a 5-ml BD Falcon polystyrene round bottom tube.

Note *Polystyrene tubes yield superior results compared to tubes composed of polyethylene or polypropylene.*

2. Add 50 μl of solution 1 and 500 μl of solution 2 to the DNA. Gently mix the contents of the tube by tapping the sides.
3. Incubate the DNA suspension at room temperature for 10–20 minutes. During this incubation period, prepare the cells as indicated in the following section (see *Preparing the Cells for Microtransfection*).

Preparing the Cells for Microtransfection

1. Remove the culture medium from the wells of the microtiter plate using the gentle vacuum aspiration of a multichannel pipettor.
2. Wash the cells by adding 200 μ l of PBS per well. Remove the PBS by aspiration.
3. To each well, add 200 μ l of standard culture medium containing 4–7% solution 3.

Note *The optimal concentration of the MBS in solution 3 should be determined by titration. A good starting point is 6%.*

The cells are now ready to receive the DNA.

Applying the DNA Suspension and Incubating the Transfection Reactions

1. Transfer 10–20 μ l of the DNA suspension into each well and incubate the microtiter plate for 3 hours at 37°C under 5% CO₂.

Notes *Incubation at 37°C under 5% CO₂ works well; however, for optimal transfection efficiency, the 3-hour incubation should be performed at 35°C under 3% CO₂.*

The optimal DNA concentration will vary, depending on the cell line and plasmid used. If higher DNA concentrations are required for optimal results, use more DNA when preparing the suspension, rather than adding more than 20 μ l of the DNA suspension to each well.

2. Remove the standard culture medium from the wells of the microtiter plate using the gentle vacuum aspiration of a multichannel pipettor.
3. Wash the cells gently with 200 μ l of PBS per well, add 200 μ l of standard growth medium into each well and incubate the microtiter plate overnight at 37°C under 5% CO₂.
4. Perform an assay to detect the gene product. For the pCMV β -gal control plasmid, perform an assay to detect the β -galactosidase gene product in order to monitor the transfection efficiency. (See *β -Galactosidase Assay for 96-Well Microtiter Plates*).

β -GALACTOSIDASE ASSAY FOR 60-MM DISHES

Harvesting the Transfected Cells

1. Gently wash the transfected cells twice with PBS.
2. Add 1.5 ml of PBS per 60-mm tissue culture dish.
3. Gently remove the cells from each tissue culture dish by scraping with a PTFE-coated scraper. Pipet the cells from each dish into a separate 1.5-ml microcentrifuge tube.
4. Pellet the cells by spinning in a microcentrifuge at $12,000 \times g$ for 5 minutes. Remove the PBS by aspiration.
5. Resuspend the cell pellet in 200 μ l of lysis buffer.[§]
6. Incubate the tubes for 5 minutes at room temperature.
7. Spin the tubes in a microcentrifuge at $12,000 \times g$ for 5 minutes.
8. Transfer the supernatant to a fresh microcentrifuge tube on ice.
9. Assay 100 μ l of cell extract for β -galactosidase activity using the protocol described in the following section.

Performing the β -Galactosidase Assay

Notes *We recommend performing the β -galactosidase assay using the same reagents supplied in the β -Galactosidase Assay Kit.*

Add fresh β -mercaptoethanol to buffer A[§] prior to the β -galactosidase assay.

1. Add 100 μ l of the cell extract to 900 μ l of buffer A.
2. Incubate for 5 minutes at 37°C.
3. Add 200 μ l of *o*-nitrophenyl- β -D-galactopyranoside (ONGP) substrate solution[§] and vortex.
4. Incubate the mixture at 37°C. Record the time and monitor the reaction until a bright yellow color develops.

Note *Incubation times for color development will vary.*

[§] See *Preparation of Media and Reagents*.

5. Terminate the color reaction by adding 500 μl of stop solution (See *Preparation of Media and Reagents*) and record the duration of the incubation (time elapsed from step 4) to use in calculating the enzyme activity.
6. Read the optical density at 420 nm (OD_{420}). (The optimum OD_{420} is between 0.6 and 0.9.)

β -GALACTOSIDASE ASSAY FOR 96-WELL MICROTITER PLATES

Note *We recommend performing the β -galactosidase assay for 96-well microtiter plates using the same reagents supplied in the β -Galactosidase Assay Kit.*

1. Remove the standard culture medium by aspiration and wash the cells once with 200 μl of PBS per well.
2. Add 25 μl of lysis buffer per well and incubate the plates for 5 minutes at room temperature.
3. Add 145 μl of buffer A per well.
4. Add 25 μl of ONGP substrate solution per well.
5. Incubate the microtiter plate at 37°C. Record the time and monitor the reaction until a bright yellow color develops (from 30 minutes to 2 hours).
6. Terminate the color reaction by adding 90 μl of stop solution (avoid creating bubbles). Record the duration of the incubation (time elapsed from step 5). Scan the microtiter plate in a microplate reader set at 414-420 nm.

β -GALACTOSIDASE HISTOCHEMICAL STAINING ASSAY

Caution *To avoid inhaling vapors, stain the cells in a fume hood.*

Notes *Perform the staining protocol within 24–72 hours after transfection.*

The quantities of the reagents given in this protocol are optimized for a 60-mm tissue culture dish.

1. Aspirate the medium from the cells.
2. Add 4 ml of freshly prepared fixing solution[§] to the tissue culture dish and incubate the dish for 10 minutes at room temperature.
3. Remove the fixing solution from the dish and gently wash the cells twice with 5 ml of PBS.
4. Add 2 ml of freshly prepared histochemical staining solution[§] to the tissue culture dish.
5. Incubate the cells between 15 minutes and overnight at 37°C in a humidified incubator.

Note *The optimal time of incubation depends on the cell type and the transfection efficiency. Observe the intensity of the blue stain in the cells under a microscope and adjust the incubation time accordingly.*

6. Remove the staining solution and wash the cells two or three times with 5 ml of PBS.
7. Add 2 ml of PBS to the tissue culture dish.

Note *For long-term storage, cover the dish with a glycerol-based mounting medium and store at 4°C.*

8. Analyze the dish under a light microscope to determine the transfection efficiency. Count the stained and unstained cells in randomly selected fields. The transfection efficiency is the percentage of stained cells in the total population.

[§] See *Preparation of Media and Reagents*

GENERAL TISSUE CULTURE TECHNIQUES: TRANSFECTION OF MAMMALIAN CELLS

Culture Inoculation

To ensure the greatest transfection efficiency, minimize trauma to the cells during incubation. For adherent cells, minimize the length of trypsinization. Remove all serum from the culture by washing with PBS, washing the culture with a medium without serum or washing the culture very quickly with a trypsin solution (serum inhibits the action of trypsin). Warm a minimal volume of trypsin to 37°C before applying to the culture. Minimal volumes for specific culture vessels are indicated in the following table:

Size and type of culture vessel	Minimum volumes
60-mm tissue culture dish	0.5 ml
100-mm tissue culture dish	1 ml
150-mm tissue culture dish	2 ml
75-cm ² tissue culture flask	1.5 ml

Proceed with trypsinization until the cell pseudopods just begin to recede toward the nucleus, a process that takes from ~30 seconds to 2 minutes, depending on the cell line. Free the cells from the dish surface by tapping the dish firmly and then immediately dilute to the desired density using medium with serum.

The best transfection efficiencies are obtained when transfected cells are grown exponentially. To allow the cells to grow as rapidly as possible, avoid contact inhibition, which occurs when the cells actually touch. Therefore, it is best to inoculate the cultures at a density that will cover no more than 15–30% of the culture vessel surface at the time of transfection.

The optimal cell density should be determined by transfecting cells plated at concentrations between 5×10^4 and 5×10^5 per 60-mm tissue culture dish.

TROUBLESHOOTING

Observation	Suggestion
The precipitate is not visible or is too fine after adding the DNA and incubating	Too much DNA was added. Check the OD ₂₆₀ of the DNA and adjust the volume used or the concentration appropriately. Transfect with different amounts of DNA to determine the optimal concentration
The precipitate is too dense, with aggregated clumps	Too little DNA was added. Check the OD ₂₆₀ of the DNA and adjust the volume used or the concentration appropriately. Transfect with different amounts of DNA to determine the optimal concentration
Unable to remove the precipitate from the plate	Ensure that the PBS does not contain Ca ²⁺ or Mg ²⁺ . It is important that the PBS wash solution contain only the cations contributed by the components listed in the recipe.

PREPARATION OF MEDIA AND REAGENTS

<p>PBS</p> <p>137 mM NaCl 2.6 mM KCl 10 mM Na₂HPO₄ 1.8 mM KH₂PO₄ Adjust the pH to 7.4 with HCl</p>	<p>Lysis Buffer</p> <p>0.25 M Tris-HCl[±] (pH 7.8) 0.5% NP40</p>
<p>Buffer A</p> <p>100 mM sodium phosphate (pH 7.5) 10 mM KCl 1 mM MgSO₄</p> <p>Note <i>Add fresh β-mercaptoethanol to buffer A prior to the assay</i></p> <p>50 mM β-mercaptoethanol</p>	
<p>ONGP Substrate Solution</p> <p>4 mg/ml in 100 mM sodium phosphate (pH 7.5)</p>	<p>Stop Solution</p> <p>1 M Na₂CO₃</p>
<p>Fixing Solution (100 ml)</p> <p>Note <i>Always use formaldehyde in a fume hood</i></p> <p>5.41 ml of 37% formaldehyde (2% final concentration) 2.00 ml of 10% glutaraldehyde (0.2% final concentration) 10.00 ml of 10× PBS (1× final concentration) Add sterile water to a final volume of 100 ml</p>	
<p>X-Gal Stock Solution</p> <p>40 mg/ml in DMSO</p>	
<p>Histochemical Reaction Mixture (250 ml)¹⁵</p> <p>0.411 g of potassium ferricyanide [K₃Fe(Cn)₆] (5 mM final concentration) 0.53 g of potassium ferrocyanide [K₄Fe(Cn)₆] (5 mM final concentration) 5 ml of 100 mM MgCl₂ (2 mM final concentration) 25 ml of 10× PBS (1× final concentration) Add X-gal just before use [final concentration of 1 mg/ml from the X-gal stock solution (see above)] Add sterile water to a final volume of 250 ml</p>	

REFERENCES

1. Chen, C. and Okayama, H. (1987) *Mol Cell Biol* 7(8):2745-52.
2. Bacchetti, S. and Graham, F. L. (1977) *Proc Natl Acad Sci U S A* 74(4):1590-4.
3. Maitland, N. J. and McDougall, J. K. (1977) *Cell* 11(1):233-41.
4. Wigler, M., Silverstein, S., Lee, L. S., Pellicer, A., Cheng, Y. *et al.* (1977) *Cell* 11(1):223-32.
5. Wigler, M., Sweet, R., Sim, G. K., Wold, B., Pellicer, A. *et al.* (1979) *Cell* 16(4):777-85.
6. Mantei, N., Boll, W. and Weissmann, C. (1979) *Nature* 281(5726):40-6.
7. Carsaro, C. M. and Pearson, M. L. (1981) *Somatic Cell Genet* 7:603-616.
8. Gorman, C., Padmanabhan, R. and Howard, B. H. (1983) *Science* 221(4610):551-3.
9. Lewis, W. H., Srinivasan, P. R., Stokoe, N. and Siminovitch, L. (1980) *Somatic Cell Genet* 6(3):333-47.
10. Luthman, H. and Magnusson, G. (1983) *Nucleic Acids Res* 11(5):1295-308.
11. Neumann, E., Schaefer-Ridder, M., Wang, Y. and Hofschneider, P. H. (1982) *Embo J* 1(7):841-5.
12. Parker, B. A. and Stark, G. R. (1979) *J Virol* 31(2):360-9.
13. Spear, D., Edwards, P., Huvar, A. and Lernhardt, W. (1993) *Strategies* 6(1):10-12.
14. Graham, F. L. and van der Eb, A. J. (1973) *Virology* 52(2):456-67.
15. Sanes, J. R., Rubenstein, J. L. and Nicolas, J. F. (1986) *Embo J* 5(12):3133-42.

MSDS INFORMATION

Material Safety Data Sheets (MSDSs) are provided online at <http://www.genomics.agilent.com>. MSDS documents are not included with product shipments.

MBS Mammalian Transfection Kit

Catalog #200388

QUICK-REFERENCE PROTOCOLS

Transfection Protocol

- ◆ Inoculate 60-mm tissue culture dishes with $1-5 \times 10^5$ exponentially growing cells/dish and grow the cells in 5 ml of appropriate medium for 24 hours
- ◆ Dilute the optimal amount of DNA (for duplicate reactions) with sterile water to a final volume of 450 μ l
- ◆ Add 50 μ l of solution 1 and 500 μ l of solution 2
- ◆ Incubate the DNA suspension at room temperature for 10–20 minutes
- ◆ Remove the culture medium from the cells by aspiration
- ◆ Wash the cells twice with PBS
- ◆ Add 5 ml of culture medium containing 4–7% solution 3
- ◆ (The cells are now ready to receive the DNA)
- ◆ Gently resuspend any precipitate in the DNA suspension with a pipettor
- ◆ Slowly add 500 μ l (1–10 μ g) of the DNA suspension to each tissue culture dish **dropwise in a circular motion**
- ◆ Swirl the tissue culture dish **once**
- ◆ Incubate for 3 hours at 37°C under 5% CO₂
- ◆ Remove the culture medium by aspiration and gently wash the tissue culture dishes three times with PBS
- ◆ Add 5 ml of standard culture medium and incubate overnight at 37°C under 5% CO₂
- ◆ Harvest the cells
- ◆ Perform an assay to detect the gene product

Microtransfection Protocol

- ◆ Inoculate a 96-well microtiter plate with $1-2 \times 10^4$ cells/well in 0.2 ml of standard culture medium and grow the cells for 24 hours
- ◆ Dilute 5–20 μ g of DNA with sterile water to a final volume of 450 μ l
- ◆ Add 50 μ l of solution 1 and 500 μ l of solution 2
- ◆ Incubate the DNA suspension at room temperature for 10–20 minutes
- ◆ Remove the culture medium from the wells of the microtiter plate using the gentle vacuum aspiration of a multichannel pipettor
- ◆ Wash the cells by adding 200 μ l of PBS/well and remove the PBS by aspiration
- ◆ Add 200 μ l of standard culture medium containing 4–7% solution 3/well
- ◆ (The cells are now ready to receive the DNA)
- ◆ Transfer 10–20 μ l (0.05–0.4 μ g) of the DNA suspension into each well
- ◆ Incubate the microtiter plate for 3 hours at 37°C under 5% CO₂
- ◆ Remove the culture medium from the wells using the gentle vacuum aspiration of a multichannel pipettor
- ◆ Wash the cells gently with 200 μ l of PBS/well
- ◆ Add 200 μ l of standard growth medium into each well and incubate overnight at 37°C under 5% CO₂
- ◆ Perform an assay to detect the gene product