



# Determination of Semivolatile Organic Compounds in Drinking Water Using the Agilent 5977A Series GC/MSD

## Application Note

Environmental

### Author

Jia-jia Wu  
Agilent Technologies Co. Ltd (China)  
Beijing,  
People's Republic of China

### Abstract

A method has been developed on the Agilent 5977A Series GC/MSD for the detection of 30 SVOCs at levels as low as 5 parts per billion (ppb), with excellent linearity ( $R^2 > 0.990$ ) and reproducibility (most relative standard deviations (RSDs)  $\leq 3.1\%$ ).

### Introduction

Semivolatile organic compounds (SVOCs) include a variety of chemicals which may have adverse health effects. Many are suspected carcinogens, and some, such as benzo(a)pyrene, are known to be carcinogenic in animals and humans. SVOCs are typically industrial plasticizers, byproducts of incomplete combustion of fossil fuels, medicines, disinfectants, pesticides, and so forth. The pollution of water can happen through contamination of the raw supply, or through leaching from sources such as landfills.

In China, regulation GB5749-2006 governs the allowed limits of SVOCs in drinking water. This application note demonstrates the ability of the 5977A Series GC/MSD to enable sensitive and reproducible detection, in 15 minutes, of 30 SVOCs covered by GB5749-2006, with RSDs that are  $\leq 3.1\%$  for most of the 30 compounds tested.



**Agilent Technologies**

## Experimental

### Standards and Reagents

Standards for 30 SVOCs and three internal standard compounds were commercially obtained (Table 1). The working calibration standards were prepared in ethyl acetate using the Agilent 7696A Sample Prep WorkBench.

Table 1. 30 SVOC Standards and Three Internal Standards (IS)

| No.  | Name                   | CAS        | No.  | Name                      | CAS        |
|------|------------------------|------------|------|---------------------------|------------|
| 1    | 1,3,5-Trichlorobenzene | 108-70-3   | 16   | Methyl parathion          | 298-00-0   |
| 2    | 1,2,4-Trichlorobenzene | 120-82-1   | 17   | Heptachlor                | 76-44-8    |
| 3    | 1,2,3-Trichlorobenzene | 87-61-6    | 18   | Malathion                 | 121-75-5   |
| 4    | Dichlorvos             | 62-73-7    | 19   | Dursban                   | 2921-88-2  |
| 5    | 2,4,6-Trichlorophenol  | 88-06-2    | 20   | Parathion                 | 56-38-2    |
| IS-1 | Acenaphthene-d10       | 15067-26-2 | 21   | Bentazone                 | 25057-89-0 |
| 6    | BHC- <i>alpha</i>      | 319-84-6   | 22   | DDE-o,p'                  | 3424-82-6  |
| 7    | Hexachlorobenzene      | 118-74-1   | 23   | DDE-p,p'                  | 72-55-9    |
| 8    | Dimethoate             | 60-51-5    | 24   | DDD-o,p'                  | 53-19-0    |
| 9    | Carbofuran             | 1563-66-2  | 25   | DDD-p,p'                  | 72-54-8    |
| 10   | Atrazine               | 1912-24-9  | 26   | DDT-o,p'                  | 789-02-6   |
| 11   | BHC- <i>beta</i>       | 319-85-7   | 27   | DDT-p,p'                  | 50-29-3    |
| 12   | Pentachlorophenol      | 87-86-5    | IS-3 | Chrysene-d12              | 1719-03-5  |
| 13   | Lindane                | 58-89-9    | 28   | Di(2-ethylhexyl)phthalate | 117-81-7   |
| 14   | Chlorothalonil         | 1897-45-6  | 29   | Benzo(a)pyrene            | 50-32-8    |
| 15   | BHC- <i>delta</i>      | 319-86-8   | 30   | Deltamethrin              | 52918-63-5 |
| IS-2 | Phenanthrene-d10       | 1517-22-2  |      |                           |            |

Table 2. Agilent 7890B GC and Agilent 5977A Series GC/MSD Conditions

#### GC run conditions

|                           |  |
|---------------------------|--|
| Analytical column         | Agilent DB-8270D UI 20 m × 180 µm, 0.36 µm (p/n 621-9723)  |
| Injection volume          | 1 µL   |
| Injection mode            | Splitless  |
| Inlet temperature         | 280 °C   |
| Liner                     | Liner, UI, splitless, single taper, no glass wool (p/n 5190-2292)                                      |
| Carrier gas               | Helium, constant flow, 1.2 mL/min  |
| Oven program              | 45 °C for 1 minute<br>50 °C/min to 100 °C<br>0.5 minutes hold<br>25 °C/min to 310 °C<br>4 minutes hold |
| Transfer line temperature | 280 °C   |

#### MS conditions

|                        |                      |
|------------------------|----------------------|
| Solvent delay          | 4 minutes            |
| Acquisition mode       | Synchronous Scan/SIM |
| Tune                   | Etune.u              |
| Gain factor            | 5.00                 |
| Source temperature     | 250 °C               |
| Quadrupole temperature | 150 °C               |

## Instruments

The study was performed on an Agilent 7890B gas chromatograph equipped with a Split/Splitless Inlet and coupled to an Agilent 5977A Series GC/MSD, using Synchronous Scan/SIM and Electron Ionization (EI) acquisition modes. Table 2 lists the instrument conditions.

## Acquisition Parameters

Table 3 shows the ions used for acquisition.

Table 3. Acquisition Parameters

| Target compounds       | RT   | Q <sub>0</sub> | Q <sub>1</sub> | Q <sub>2</sub> | Q <sub>3</sub> | Target compounds          | RT    | Q <sub>0</sub> | Q <sub>1</sub> | Q <sub>2</sub> | Q <sub>3</sub> |
|------------------------|------|----------------|----------------|----------------|----------------|---------------------------|-------|----------------|----------------|----------------|----------------|
| 1,3,5-Trichlorobenzene | 4.75 | 180            | 182            | 184            | 145            | Methyl parathion          | 8.67  | 109            | 125            | 263            | 79             |
| 1,2,4-Trichlorobenzene | 5.05 | 180            | 182            | 184            | 145            | Heptachlor                | 8.82  | 100            | 272            | 274            | 270            |
| 1,2,3-Trichlorobenzene | 5.26 | 180            | 182            | 184            | 145            | Malathion                 | 8.9   | 173            | 127            | 125            | 93             |
| Dichlorvos             | 5.34 | 109            | 185            | 79             | 187            | Dursban                   | 8.99  | 197            | 199            | 314            | 97             |
| 2,4,6-Trichlorophenol  | 6.06 | 196            | 198            | 97             | 132            | Parathion                 | 9.06  | 109            | 97             | 291            | 125            |
| BHC- <i>alpha</i>      | 7.88 | 219            | 183            | 181            | 217            | Bentazone                 | 9.15  | 198            | 119            | 161            | 92             |
| Hexachlorobenzene      | 7.92 | 284            | 286            | 282            | 288            | DDE-o,p'                  | 9.59  | 246            | 248            | 318            | 316            |
| Dimethoate             | 7.95 | 87             | 93             | 125            | 143            | DDE-p,p'                  | 9.84  | 246            | 248            | 318            | 316            |
| Carbofuran             | 7.96 | 164            | 149            | 131            | 123            | DDD-o,p'                  | 9.91  | 235            | 237            | 165            | 212            |
| Atrazine               | 8.01 | 200            | 215            | 58             | 173            | DDD-p,p'                  | 10.17 | 235            | 237            | 165            | 212            |
| BHC- <i>beta</i>       | 8.07 | 219            | 181            | 183            | 109            | DDT-o,p'                  | 10.21 | 235            | 237            | 165            | 212            |
| Pentachlorophenol      | 8.11 | 266            | 264            | 268            | 165            | DDT-p,p'                  | 10.47 | 235            | 237            | 165            | 212            |
| Lindane                | 8.17 | 219            | 183            | 181            | 111            | Di(2-ethylhexyl)phthalate | 10.94 | 149            | 167            | 57             | 70             |
| Chlorothalonil         | 8.29 | 266            | 264            | 268            | 270            | Benzo(a)pyrene            | 12.61 | 252            | 253            | 125            | 126            |
| BHC- <i>delta</i>      | 8.39 | 219            | 183            | 181            | 217            | Deltamethrin              | 13.29 | 181            | 253            | 251            | 255            |

RT – retention time in minutes

Q<sub>0</sub> – quantifier ion

Q<sub>1</sub>, Q<sub>2</sub>, Q<sub>3</sub> – qualifier ions

Table 4. Calibration Coefficients for the 30 SVOCs

| Target compound        | R <sup>2</sup> | Target compound           | R <sup>2</sup> |
|------------------------|----------------|---------------------------|----------------|
| 1,3,5-Trichlorobenzene | 0.999          | Methyl parathion          | 0.994          |
| 1,2,4-Trichlorobenzene | 0.999          | Heptachlor                | 0.997          |
| 1,2,3-Trichlorobenzene | 0.999          | Malathion                 | 0.993          |
| Dichlorvos             | 0.998          | Dursban                   | 0.991          |
| 2,4,6-Trichlorophenol  | 0.995          | Parathion                 | 0.993          |
| BHC- <i>alpha</i>      | 0.999          | Bentazone**               | 0.991          |
| Hexachlorobenzene      | 0.999          | DDE-o,p'                  | 0.999          |
| Dimethoate             | 0.995          | DDE-p,p'                  | 0.999          |
| Carbofuran             | 0.993          | DDD-o,p'                  | 0.999          |
| Atrazine               | 0.994          | DDD-p,p'                  | 0.996          |
| BHC- <i>beta</i>       | 0.999          | DDT-o,p'                  | 0.998          |
| Pentachlorophenol*     | 0.993          | DDT-p,p'                  | 0.993          |
| BHC- <i>gamma</i>      | 0.999          | Di(2-Ethylhexyl)phthalate | 0.997          |
| Chlorothalonil         | 0.997          | Benzo(a)pyrene            | 0.998          |
| BHC- <i>delta</i>      | 0.999          | Deltamethrin              | 0.999          |

\*pentachlorophenol: 20–200 ppb

\*\* bentazone: 50–1,000 ppb

## Results and Discussion

### Linearity

Calibration curves were constructed from 5 to 200 ng/mL ppb for most target compounds, using three internal standards to ensure coverage of all compounds across the entire retention time range. Figure 1 shows the curves for four of the compounds, which had calibration coefficient ( $R^2$ ) values  $\leq 0.997$ . Pentachlorophenol curves were constructed from 20–200 ng/mL, and bentazone calibration curves were constructed from 50–1,000 ng/mL. The  $R^2$  values for pentachlorophenol and bentazone were greater than 0.990.




Figure 1. Example calibration curves for four of the SVOCs, from 5 to 200 ng/mL.

## Reproducibility

Table 5 illustrates the excellent reproducibility of quantitation obtainable with the 5977A Series GC/MSD. Using the Etune setting provided the lowest RSDs across all 30 SVOC compounds, most of which were  $\leq 3.1\%$  across eight injections, and seven were less than 1%.

Table 5. Reproducibility (RSD) for a 40 ppb Standards Mix Sample\*

| Target compound        | Average concentration (ppb) | RSD (%) | Target compound           | Average concentration (ppb) | RSD (%) |
|------------------------|-----------------------------|---------|---------------------------|-----------------------------|---------|
| 1,3,5-Trichlorobenzene | 40.5                        | 0.7     | Methyl parathion          | 39.9                        | 1.5     |
| 1,2,4-Trichlorobenzene | 40.5                        | 0.5     | Heptachlor                | 37.7                        | 0.8     |
| 1,2,3-Trichlorobenzene | 40.5                        | 0.7     | Malathion                 | 36.6                        | 1.2     |
| Dichlorvos             | 38.8                        | 4.6     | Dursban                   | 36.4                        | 1.4     |
| 2,4,6-Trichlorophenol  | 38.9                        | 1.7     | Parathion                 | 35.0                        | 2.2     |
| BHC- <i>alpha</i>      | 38.3                        | 1.1     | Bentazone                 | 173.0                       | 2.7     |
| Hexachlorobenzene      | 40.6                        | 1.0     | DDE-o,p'                  | 37.7                        | 0.6     |
| Dimethoate             | 35.5                        | 3.7     | DDE-p,p'                  | 37.4                        | 1.0     |
| Carbofuran             | 32.8                        | 3.5     | DDD-o,p'                  | 36.4                        | 0.9     |
| Atrazine               | 35.9                        | 1.0     | DDD-p,p'                  | 35.0                        | 1.4     |
| BHC- <i>beta</i>       | 38.1                        | 0.7     | DDT-o,p'                  | 34.4                        | 1.2     |
| Pentachlorophenol      | 41.2                        | 3.1     | DDT-p,p'                  | 37.1                        | 4.4     |
| BHC- <i>gamma</i>      | 37.8                        | 1.2     | Di(2-ethylhexyl)phthalate | 37.5                        | 3.1     |
| Chlorothalonil         | 34.9                        | 1.5     | Benzo(a)pyrene            | 33.5                        | 2.3     |
| BHC- <i>delta</i>      | 38.3                        | 1.5     | Deltamethrin              | 37.5                        | 2.0     |

\* 8 consecutive injections were used to calculate the RSDs. The concentration of bentazone was 200 ppb.

## The Etune Setting Provides Better Sensitivity

The Extractor Ion Source and tuning protocol (Etune) have increased the sensitivity of the 5977A Series GC/MSD. For example, using Etune versus Atune with a 40 ppb SVOC standards mix provides at least two-fold higher peak heights for every compound in the mix (Figure 2).




Figure 2. Extracted ion chromatogram for a 40 ppb standards mix solution using the Atune (red) or Etune (blue) tuning protocol, with one region expanded to illustrate the difference in peak heights.

## **Conclusion**

The Agilent 5977A Series GC/MSD provides a stable platform for the sensitive, accurate, and reproducible analysis of SVOCs in drinking water, to meet the Chinese regulation. Use of the Etune tuning protocol provides the highest reproducibility and sensitivity. Using automation such as the Agilent 7696A Sample Prep WorkBench to prepare standards also results in reproducible calibration coefficient ( $R^2$ ) values.

## **For More Information**

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2013  
Printed in the USA  
April 3, 2013  
5991-2127EN



**Agilent Technologies**