

**Dusíkatá
rozpouštědla
(Amidy)**

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Dimethyl Formamide	✓	✓	✗	✗	✗	✓	✓	✓
Diethylacetamide	✓	✓	-	-	✗	✓	-	✓
Triethanolamine	✓	✓	-	-	✓	✓	-	-
Aniline	-	✓	-	-	✗	✓	-	-
Pyridine	✓	✓	✓	✗	✗	✓	✗	✓
Acetonitrile	✓	✓	✓	✓	✗	✓	✓	✓

Uhlovodíky

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Hexane, Xylene	✓	✓	✓	✗	✓	✓	✗	✓
Toluene, Benzene	✓	✓	✓	✗	✓	✓	✗	✓
Kerosene, Gasoline	✓	✓	✓	✓	✓	✓	✓	-
Tetralin, Decalin	-	✓	✓	-	✓	✓	-	-

**Halogenované
uhlovodíky**

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Methylene Chloride	✓	✓	✓	✗	✗	✓	✓	✓
Chloroform	✓	✓	✓	✗	✗	✓	✓	✓
Trichloroethylene	✓	✓	✓	✗	✓	✓	✓	✓
Monochlorobenzene, Freon	✓	✓	✓	✓	✓	✓	✓	✓
Carbon Tetrachloride	✓	✓	✓	✗	✓	✓	✓	✓

Ketony

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Acetone, Cyclohexanone	✓	✓	✗	✗	✗	✓	✓	✓
Methyl Ethyl Ketone	✓	✓	✓	✗	✓	✓	✓	✓
Isopropylacetone	✓	✓	✗	✗	✓	✓	-	✓
Methyl Isobutyl Ketone	-	✓	✓	✗	-	✓	✓	✓

Estery

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Ethyl & Methyl Acetate	✓	✓	✓	✗	✗	✓	✓	✓
Amyl & Butyl Acetate	✓	✓	✗	✗	✓	✓	✓	✓
Propyl Acetate	✓	✓	✗	✗	✓	✓	✓	-
Propylene Glycol Acetate	-	✓	-	✗	✗	✓	✓	-
2-Ethoxyethyl Acetate	-	✓	-	✗	✓	✓	-	-
Methyl Cellosolve Acetate	-	✓	-	✗	✗	✓	✓	✓
Benzyl Benzoate	✓	✓	-	✗	✓	✓	-	-
Isopropyl Myristate	✓	✓	-	✗	✓	✓	-	-
Tricresyl Phosphate	-	✓	-	✗	✓	✓	-	-

**Oxidy –
Ethery**

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Ethyl Ether	✓	✓	✓	✓	✓	✓	✓	-
Dioxane	✓	✓	✓	✗	✗	✓	-	✓
Tetrahydrofuran	✓	✓	✓	✗	✗	✓	-	✓
Triethanolamine	✓	✓	-	-	✓	✓	-	-
Dimethylsulfoxide (DMSO)	✓	✓	✗	✗	✗	✓	✗	✓
Isopropyl Ether	-	✓	✓	✓	✓	✓	✓	-

Ostatní

	Nylon	PTFE	PVDF	PES	CA	RC	PP	GMF
Phenol, Aqueous, 10%	-	✓	✓	✗	✗	✗	✓	✓
Formaldehyde Solution, 30%	✓	✓	✓	✓	✓	✓	✓	✓
Hydrogen Peroxide, 30%	✓	✓	-	-	✓	✓	-	-
Silicone Oil & Mineral Oil	-	✓	✓	✓	✓	✓	✓	✓