

## MS-PA LED Digital Magnetic Stirrer

### Features

- Digital speed control within the range of 0-1500rpm
- Max. stirring quantity of H<sub>2</sub>O at 3L
- LED display shows speed
- Nylon+GF housing provides a good performance of chemical resistance



## MS-PB Magnetic Stirrer

### Features

- A wide range of 0-1500rpm
- Max. stirring quantity of H<sub>2</sub>O at 3L
- Nylon+GF housing provides a good performance of chemical resistance



Chemical resistance



Strong magnetism



LED display



### Specifications

Work plate dimension	
Work plate material	
Motor type	
Motor input power	
Motor output power	
Power	
Voltage	
Stirring positions	
Max. stirring quantity[H <sub>2</sub> O]	
Max. magnetic bar[length]	
Speed range	
Speed display	
Protection class	
Dimension[WxDxH]	
Weight	
Permissible ambient temperature and humidity	

### MS-PA

φ135mm(5 inch)  
Nylon+GF  
Brushless DC motor

5W

3W

15W

100-120/200-240V,50/60Hz

1

3L

50mm

0-1500rpm

LED

IP42

150×260×80mm

1.8kg

10-40°C, 80%RH

### MS-PB

φ135mm(5 inch)  
Nylon+GF  
Brushless DC motor

5W

3W

10W

100-120/200-240V,50/60Hz

1

3L

50mm

0-1500rpm

Scale

IP42

150×260×80 mm

1.8kg

10-40°C, 80%RH