




MANUFACTURER OF LABORATORY AND UTILITY GLASSWARE

Petri dishes ANUMBRA® and watch glass


FGH plus, s.r.o., 8. května 33, 787 01 Šumperk, CZ

PETRI DISH
made of soda lime glass



DIMENSIONS IN MM

SET MARKING	LID		BOTTOM		CURVATURE MAX.	GLASS THICKNESS t
	d 1	h 1	d 2	h 2		
40 x 12	34,5 + 1,5 - 0,0	8,4 + 1,6 - 2,4	33,0 + 1,0 - 1,0	12,0 + 1,0 - 1,0	1,2	1,9 + 0,2 - 0,2
50 x 12	45,0 + 1,5 - 0,0	8,4 + 1,6 - 2,4	44,0 + 1,0 - 1,0	12,0 + 1,0 - 1,0	1,2	1,9 + 0,2 - 0,2
60 x 12	55,5 + 1,5 - 0,0	8,4 + 1,6 - 2,4	54,0 + 1,0 - 1,0	12,0 + 1,2 - 1,2	1,2	1,9 + 0,2 - 0,2
60 x 15	55,5 + 1,8 - 0,0	12,4 + 1,6 - 2,4	54,0 + 1,0 - 1,0	15,0 + 1,3 - 1,3	1,2	1,9 + 0,2 - 0,2
70 x 15	65,5 + 1,8 - 0,0	12,4 + 1,6 - 2,4	64,0 + 1,0 - 1,0	15,0 + 1,3 - 1,3	1,2	1,9 + 0,2 - 0,2
80 x 15	75,5 + 1,8 - 0,0	12,4 + 1,6 - 2,4	74,0 + 1,0 - 1,0	15,0 + 1,3 - 1,3	1,2	1,9 + 0,2 - 0,2
90 x 15	85,5 + 2,0 - 0,0	12,4 + 1,6 - 2,4	84,0 + 1,0 - 1,0	15,0 + 1,3 - 1,3	1,3	1,9 + 0,2 - 0,2
100 x 10	95,5 + 2,0 - 0,0	6,4 + 1,6 - 2,4	94,0 + 1,0 - 1,0	10,0 + 1,2 - 1,2	1,4	1,9 + 0,2 - 0,2
100 x 15	95,5 + 2,0 - 0,0	12,4 + 1,6 - 2,4	94,0 + 1,0 - 1,0	15,0 + 1,3 - 1,3	1,4	1,9 + 0,2 - 0,2
100 x 20	95,5 + 2,0 - 0,0	15,4 + 1,6 - 2,4	94,0 + 1,0 - 1,0	20,0 + 1,5 - 1,5	1,4	1,9 + 0,2 - 0,2
120 x 20	115,5 + 2,0 - 0,0	15,4 + 1,6 - 2,4	114,0 + 1,0 - 1,0	20,0 + 1,5 - 1,5	1,4	1,9 + 0,2 - 0,2
150 x 25	143,0 + 2,0 - 0,5	18,4 + 1,9 - 2,7	141,0 + 1,0 - 1,5	25,0 + 1,8 - 1,8	1,5	1,9 + 0,2 - 0,2
180 x 30	172,5 + 2,0 - 0,5	24,4 + 1,9 - 2,7	171,0 + 1,0 - 1,5	30,0 + 2,0 - 2,0	1,5	1,9 + 0,2 - 0,2
200 x 30	194,5 + 2,0 - 1,0	24,4 + 2,1 - 2,9	191,5 + 1,5 - 1,5	30,0 + 2,0 - 2,0	1,6	1,9 + 0,3 - 0,3



ANUMBRA