

## Additional peristaltic pumps

### ADU (art. 902362)

#### UNIVERSAL PERISTALTIC PUMP FOR LIQUID ADDITIVES

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



#### Machines

GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA	GW3060S
GW3060SA	GW3060SC	GW3060SCA	GW4090	GW4090A
GW4090BXC	GW4090BXCA	GW4090BXP	GW4090C	GW4090CA
GW4090P	GW4090PA			

### FLUX4060 (art. 902900)

#### DOSING CONTROL FOR ADDITIONAL PERISTALTIC PUMPS

- allows double control for proper dosage of detergent and liquid neutralizer



#### Machines

GW4260-00-000	GW4260-S0-000	GW4260-00-000	GW4260-S0-000	GW4260-S0-060
GW4260-SC-000	GW1260-00-000	GW1260-S0-000	GW1260-SC-000	

### P12145 (art. 903766)

#### PERISTALTIC PUMP FOR LIQUID DETERGENT

- for dosing liquid alkaline detergent
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



#### Machines

GW1245-S0-100
---------------

**P14060** (art. 902904)

**PERISTALTIC PUMP FOR LIQUID DETERGENT**

- for dosing liquid alkaline detergent
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000	

**P14190** (art. 903418)

**PERISTALTIC PUMP FOR LIQUID DETERGENT**

- for dosing liquid alkaline detergent
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW4290-00-000	GW4290-S0-000	GW4290-S0-000		
---------------	---------------	---------------	--	--

**P34090** (art. 904832)

**PERISTALTIC PUMP FOR CAUSTIC SODA**

- for dosing caustic soda solution
- with stainless steel drawing pipe

**Machines**

GW4090	GW4090A	GW4090BXC	GW4090BXCA	GW4090BXP
GW4090C	GW4090CA	GW4090P	GW4090PA	

**P34260** (art. 905436)

**PERISTALTIC PUMP FOR CAUSTIC SODA**

- for dosing caustic soda solution
- with improved PVC drawing pipe for high chemical resistance
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000	
---------------	---------------	---------------	---------------	--

**P34290** (art. 905438)

**PERISTALTIC PUMP FOR CAUSTIC SODA**

- for dosing caustic soda solution
- with improved PVC drawing pipe for high chemical resistance
- suitable for 2 liters, 5 liters cans



**Machines**

GW4290-00-000      GW4290-S0-000

**P36010** (art. 903531)

**PERISTALTIC PUMP FOR LIQUID ADDITIVES**

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW6010T	GW6010TM	GW7010	GW7010M	GW7010T
GW7010TM	GW7015	GW7015T		

**P36290** (art. 904769)

**PERISTALTIC PUMP FOR LIQUID ADDITIVES**

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans

**Machines**

GW6290DS      GW6290DSG

**P44090** (art. 904831)

---

**PERISTALTIC PUMP FOR LIQUID ANTIFOAM**

- for dosing liquid antifoam
- with stainless steel drawing pipe

**Machines**

GW4090	GW4090A	GW4090BXC	GW4090BXCA	GW4090BXP
GW4090C	GW4090CA	GW4090P	GW4090PA	

**P44260** (art. 905440)

---

**PERISTALTIC PUMP FOR LIQUID ANTIFOAM**

- for dosing liquid antifoam
- with improved PVC drawing pipe for high chemical resistance
- suitable for 2 liters, 5 liters cans



**Machines**

GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000
---------------	---------------	---------------	---------------

**P44290** (art. 905441)

---

**PERISTALTIC PUMP FOR LIQUID ANTIFOAM**

- for dosing liquid antifoam
- with improved PVC drawing pipe for high chemical resistance
- suitable for 2 liters, 5 liters cans



**Machines**

GW4290-00-000	GW4290-S0-000
---------------	---------------

**P46010** (art. 903530)

**PERISTALTIC PUMP FOR LIQUID ANTIFOAM**

- for dosing liquid antifoam
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW6010T	GW6010TM	GW7010	GW7010M	GW7010T
GW7010TM	GW7015	GW7015T		

**P46010L** (art. 905017)

**PERISTALTIC PUMP FOR DOSING LIQUID ADDITIVES SOLUTION TO AVOID FOAM**

- Peristaltic pump for dosing liquid additives solution to avoid foam
- With stainless steel drain
- Suitable for 2 liter, 5 liter and 10 liter tanks
- Specific for GW6010L series machines

**Machines**

GW6010L	GW6010LX	GW6010ML	GW6010MLX
---------	----------	----------	-----------

**P46290** (art. 904768)

**PERISTALTIC PUMP FOR LIQUID ADDITIVES**

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans

**Machines**

GW6290DS	GW6290DSG
----------	-----------

**P56090B** (art. 903609)

---

**PERISTALTIC PUMP FOR LIQUID ADDITIVES**

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans

**Machines**

GW6290DS                      GW6290DSG

**P57015** (art. 904266)

---

**PERISTALTIC PUMP FOR LIQUID ADDITIVES**

- for dosing liquid additives (ex: alkaline detergent, caustic solution or antifoam)
- with stainless steel drawing pipe
- suitable for 2 liters, 5 liters, 10 liters cans



**Machines**

GW7015                      GW7015T

**TANKBSK** (art. 902591)

---

**UPPER LEVEL CAN HOLDER BASKET**

- capacity: 2 x 5 liters jerrycans
- necessary in case of additional pumps WD-AD5HT or WD-AD6HT
- expands detergent cabinet capacity from 2 to 4 x 5 liters jerrycans



**Machines**

GW4290-00-000              GW4290-S0-000

**WD-FLUX** (art. 901703)

**FLOW METER FOR ADDITIONAL PERISTALTIC PUMPS**

- allows redundant control for proper dosage of detergent and liquid neutralizer



**WD-FLUX6010** (art. 903459)

**DETERGENT DOSAGE CONTROL FOR ADDITIONAL PERISTALTIC PUMPS**

- allows double control for proper dosage of detergent and liquid neutralizer



**Machines**

GW6010L	GW6010LX	GW6010ML	GW6010MLX	GW6010T
GW6010TM	GW7010	GW7010M	GW7010T	GW7010TM
GW7015	GW7015T			

**WD-FLUX6290** (art. 904770)

**DETERGENT DOSAGE CONTROL FOR ADDITIONAL PERISTALTIC PUMPS**

- allows double control for proper dosage of detergent and liquid neutralizer

**Machines**

GW6290DS

**WD-LS3060** (art. 902084)

**DETERGENT LEVEL SENSOR**

- with adjustable sensor for cans of different heights
- suitable for all models from 45 - 60 cm



**Machines**

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA
GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000

**WD-LS4190** (art. 903167)

**DETERGENT LEVEL SENSOR**

- with adjustable sensor for cans of different heights
- suitable for 90 cm models



**Machines**

GW4290-00-000      GW4290-S0-000      GW4290-SC-000

**WD-LS5090** (art. 901710)

**DETERGENT LEVEL SENSOR**

- with adjustable sensor for cans of different heights
- suitable for all models 90 cm and large dimensions



**Machines**

GW4090              GW4090A              GW4090BXC              GW4090BXCA              GW4090BXP  
 GW4090C              GW4090CA              GW4090P              GW4090PA

**WD-LS6010** (art. 903525)

**DETERGENT LEVEL SENSOR**

- with adjustable sensor for cans of different heights
- suitable for 6010/7010/6290 series



**Machines**

GW6010L              GW6010LX              GW6010ML              GW6010T              GW6010TM  
 GW6290DS              GW7010              GW7010M              GW7010T              GW7010TM  
 GW7015              GW7015T