

## Booster pumps for demineralized water

### PAD (art. 900644)

#### BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER

- allows machine water supply with demineralized water from a non-pressure tank positioned at a minimum height from the floor of 85 cm



#### Machines

GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA	GW3060P
GW3060PA	GW3060S	GW3060SA	GW3060SC	GW3060SCA
GW4090	GW4090A	GW4090BXC	GW4090BXCA	GW4090BXP
GW4090C	GW4090CA	GW4090P	GW4090PA	

### PAD1 (art. 901394)

#### BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER

- allows machine water supply with demineralized water drawn from a non-pressure tank positioned at a minimum height from the floor of 85 cm



#### Machines

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000
GW4290-00-000	GW4290-S0-000			

### PAD2 (art. 900989)

#### BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER

- used on machines for supplying demineralized water from a tank placed on the floor
- could require PAD2R kit



#### Machines

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA

GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4090	GW4090A	GW4090BXC	GW4090BXCA
GW4090BXP	GW4090C	GW4090CA	GW4090P	GW4090PA
GW4260-00-000	GW4260-S0-000	GW4260-SC-000	GW4290-00-000	GW4290-SC-000
GW6290DS				

**PAD2-60 (art. 901789)**

**BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER**

- allows machine with demineralized water supply from a tank which must be positioned on the floor
- could require PAD2R kit



**Machines**

GW4260-S0-060	GW2145	GW2145A	GW3060BX	GW3060BXA
GW3060BXC	GW3060BXCA	GW3060P	GW3060PA	GW3060S
GW3060SA	GW3060SC	GW3060SCA	GW4090	GW4090A
GW4090BXC	GW4090BXCA	GW4090BXP	GW4090C	GW4090CA
GW4090P	GW4090PA	GW6290DS		

**PAD26010 (art. 903840)**

**BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER**

- used on machines for supplying demineralized water from a tank placed on the floor



**Machines**

GW6010L	GW6010LX	GW6010ML	GW6010MLX	GW6010T
GW6010TM	GW7010	GW7010M	GW7010T	GW7010TM
GW7015	GW7015T			

**PAD26010-60** (art. 904517)

**BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER - 60Hz**

- used on machines for supplying demineralized water from a tank placed on the floor



**PAD2R** (art. 903001)

**PAD2 RELAY KIT**

- controls PAD2 pump on basic machines



**Machines**

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW4260-00-000	GW4260-SC-000	GW4290-00-000	GW4290-SC-000

**PAD2X** (art. 903330)

**ST. STEEL BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER**

- used on machines for supplying demineralized water from a tank placed on the floor
- pump head made of stainless steel
- self priming
- could require PAD2R kit



**Machines**

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA
GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4090	GW4090A	GW4090BXC	GW4090BXCA
GW4090BXP	GW4090C	GW4090CA	GW4090P	GW4090PA
GW4260-00-000	GW4260-S0-000	GW4260-SC-000	GW4290-00-000	GW4290-S0-000
GW4290-SC-000	GW6290DS			

**PAD2X-60** (art. 903332)

**BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER - 60 HZ**

- allows machine with demineralized water supply from a tank which must be positioned on the floor
- self priming
- could require PAD2R kit



**Machines**

GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA
GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4090	GW4090A	GW4090BXC	GW4090BXCA
GW4090BXP	GW4090C	GW4090CA	GW4090P	GW4090PA
GW4260-S0-060	GW6290DS			

**PAD2X6010** (art. 903841)

**ST. STEEL BOOSTER PUMP FOR NON-PRESSURE DEMINERALIZED WATER**

- used on machines for suppling demineralized water from a tank placed on the floor
- self priming
- pump head made of stainless steel



**Machines**

GW6010L	GW6010LX	GW6010ML	GW6010MLX	GW6010T
GW6010TM	GW7010	GW7010M	GW7010T	GW7010TM
GW7015	GW7015T			