

Agilent

SureScan Dx Microarray Scanner

Compact, Sensitive, and Accurate

The Measure of Confidence

CE Marked for Invitro Diagnostic Use in Europe*



Agilent's CE-IVD Microarray Scanner is Ideal for Clinical Labs and Hospitals

The SureScan Dx Microarray Scanner is the foundation of our complete microarray solution and represents the latest innovation in scanner technology from Agilent.

Exceptional limits of detection provide the ability to obtain as much biological information as possible, with sensitivity and resolution — from a single data point or a single experiment.

Confidence in Results

- The precision and reproducibility study demonstrated that Gene Expression and CGH microarray probes were highly reproducible when scanned across multiple microarrays, scanners, operators and days.
- Exceptional limit of detection allows accurate measurement of very low signals.
- Precision-engineered optics enable high feature resolution.
- Built-in ozone protection minimizes signal degradation.

Streamlined Workflow

- Continuous slide loading capability to remove batch-loading restrictions
- Compact footprint optimizes bench space utilization

Learn how to gain more confidence in your results today!

www.agilent.com/genomics/SureScanDx

Reproducible Data

Gene Expression: Background subtracted signals, for both the red and green channels, had a 1-sided 95% CI of Coefficient of Variation of <15%: 7.92% in the green channel and 8.92% in the red channel.

CGH: The log ratio of CGH probes had a 1-sided 95% CI of standard deviation of <0.06 at 0.03.



Agilent Technologies

Ordering Information

Product	Product Number
SureScan Microarray Dx Scanner Bundle (For sale in registered European Union countries)	G5761AA
SureScan Microarray Scanner Bundle (Research Use Only)	G4900DA

Specifications

Feature	Description
Dynamic Range	>104 (16-bit data format), >105 (20-bit data format), 106 (with XDR scanning)
Resolution	2, 3, 5 and 10 microns
Dynamic Auto-Focus	Continually adjusts scanner's focus, keeping in focus at all times
Autoloader	24-slide cassette allows for hands-off operation with continuous access for limitless scan runs
Integrated Barcode Reader	Reads code 128, Code 39, Code 93, and CODABAR
Compatible Dyes	Cyanine 3 and Cyanine 5, and Alexa 647, 555, and 660
Laser Information	Green solid-state laser, 532 nm; Red solid-state laser, 640 nm Power: 20 mW at 532 nm and 640 nm both controlled to 13 mW
Scan Window Maximum	71 mm x 21.6 mm
PMT Adjustment	Automatic PMT gain calibration before each run allows adjustment of signal levels from 100% (default) to 1%
Detection Limit	0.01 chromophores per square micron
Pixel Placement Error	<1 pixel at 5-micron resolution
Uniformity	5% CV global nonuniformity, average local nonuniformity is typically 1% based upon 100-micron features
Scan Time	2-color simultaneous data acquisition in 16 minutes per for 3-micron scans and 24 minutes for 2-micron scans (scan region of 61 mm x 21.6 mm)
Data Workstation and Operation System	PC based with Windows 7 Professional – 64 bit
Weight Scanner	125 lbs (56.8 kg)
Power Input	100 – 240 Vac, 50 – 60 Hz, 250 VA max.
Fuses	Two power supply fuses: T4A, 250 VAC
Humidity Operating	15% to 95% RH at 30°C
Altitude	Operating maximum: 2,300 m (7,500 ft)
Operating Temperature Range	15° to 30°C
Laser Product Classification	Class 1

Intended Use: The SureScan Dx Microarray Scanner system consisting of SureScan Dx Microarray Scanner with autoloader, and Agilent Microarray Scan Control software is intended to measure fluorescence signals of labeled DNA and RNA target hybridized to microarrays.

Indications for Use: The SureScan Dx Microarray Scanner system is indicated for use in a clinical laboratory environment when measuring fluorescence signals of labeled DNA and RNA target hybridized to microarrays used as part of a validated diagnostic assay.

* The SureScan Dx Microarray Scanner is CE- marked and is currently marketed in Belgium, Denmark, Finland, France, Germany, Italy, Netherlands, Spain, Sweden, and the United Kingdom.

Learn more at:

www.agilent.com/genomics/surescan

For In Vitro Diagnostic Use. The SureScan Dx Microarray Scanner system is approved for In Vitro Diagnostic Use within the European Union

© Agilent Technologies, Inc. 2013
Published in USA, February 1, 2013
5991-1928EN

The Measure of Confidence



Agilent Technologies